


# MAGO 4S

Multiple solutions  
Infinite possibilities



- Automation of ELISA, IFA and agglutination tests processing simultaneously in one batch
- 120 sample positions
- Performs up to 12 ELISA test (on 4 microplates) or up to 16 IFA test
- Levey-Jennings and Westgards rules program on board

 **Delta Biologicals srl**  
Via Nicaragua 12/14 - 00071  
Pomezia (Roma) - Italy  
+39 06 91190232- 6+39 06 9105244



MAGO4S27102017EN

Specifications are subject to change without prior notice.  
Pictures used only for representation purposes- subject to change without prior notice.

# MAGO 4S



Food Intolerance  
& Allergy

Tumor Marker Panel

TORCH/EBV Panel

**Multiple solutions,  
infinite possibilities.**

**Autoimmune Panel**

- Celiac Disease
- Rheumatic Disease
- Vascular Diseases

**Hormones Panel**

- Fertility
- Steroids
- Thyroid

Hepatitis Panel



Separate  
built/in  
compartment  
for compressor  
and hydraulic  
system.



An independent  
incubation  
modules  
for each plate.



Automated  
well-washing  
performed  
at customized  
volumes  
and cycles.



Software  
direct placement  
of reagents  
and the amount  
of each reagent  
needed.

## Functional characteristics

- Ability to perform multiple assays per sample
- Liquid sensing mechanism reliability
- Pre-assay plate and reagent volume check
- Erba reagents in ready to use containers
- Serial two-fold dilutions performed on-board
- User friendly Windows 10 based software
- Integrated ELISA reader, halogen lamp, (interference filters 405, 450, 492, 620 (two free positions))

## Technical data

Dimensions: 120 x 70 x 60 cm  
Weight: 80 kg  
Power Supply voltage: 115-230 Volts

## Computer Characteristics

15,6 flat colour monitor  
Pentium® processor  
500 GB hard disk  
4GB RAM  
USB Port, serial port, ethernet port  
Windows 10 software platform