



Applications for Clinical Microbiology Testing
Applications pour la
microbiologie clinique

2019



BioMaxima S.A. is an European Union-based manufacturer of a wide range of in vitro diagnostic (IVD) products. The company was established in 2002 in Poland and it was listed on the Warsaw Stock Exchange in 2010.

BioMaxima products are CE marked and meet the requirements of EC directive 98/79/EC. Our Quality System is certified to the ISO 9001:2015 and ISO 13485:2016 international standards and quality control is carried out in accordance with ISO 11133:2014.

BioMaxima is focused on continuous growth and development, working internationally through a network of authorised distributors. Wherever BioMaxima products reach, our services and product support follow, meeting the highest standards. We supply international markets including countries of Europe, Asia, Africa and South America.

BioMaxima S.A. est un fabricant européen d'une large gamme de dispositifs de diagnostic in vitro (IVD). La société a été créée en 2002 en Pologne et a été cotée à la Bourse de Varsovie en 2010.

Les produits BioMaxima sont marqués CE et répondent aux exigences de la directive CE 98/79 / CE. Notre système de qualité est certifié selon les normes internationales ISO 9001: 2015 et ISO 13485: 2016, et le contrôle de la qualité est effectué conformément à la norme ISO 11133: 2014.

BioMaxima se concentre sur le développement continu, travaillant à l'international à travers un réseau de distributeurs agréés. Quel que soit la destination de nos produits nous assurons leur entretien et le service après-vente aux standards les plus exigeants. Nous fournissons des marchés internationaux comprenant des pays d'Europe, d'Asie, d'Afrique et d'Amérique du Sud.

Applications for Pharmaceutical and Cosmetics Microbiology Testing
Applications pour analyses de microbiologie pharmaceutique et cosmétique

Applications for Food Microbiology Testing
Applications pour analyses de microbiologie alimentaire

Applications for Water Microbiology Testing
Applications pour analyses de microbiologie de l'eau

Applications for Clinical Microbiology
Applications pour la microbiologie clinique

Applications for Molecular Biology
Applications pour la biologie moléculaire

IVD Products
Dispositifs médicaux de diagnostic in vitro

Product	Page
1. Prepared culture media/milieux de culture prepares	4
1.1. Prepared culture media on petri dishes/Milieux de culture prepares petri	4
1.1.1. Petri dishes Ø 90 mm	4
1.1.2. Petri dishes Ø 140 mm	7
1.1.3. Two-sections Petri dishes Ø 90 mm	7
1.1.4. Three-sections Petri dishes Ø 90 mm	9
1.1.5. Four-sections Petri dishes Ø 90 mm	9
1.1.6. Control of surfaces Petri dishes Ø 55 mm	9
1.2. Prepared culture media in the bottles / Milieux de culture prepares en bouteilles	11
1.3. Prepared culture media in tubes/Milieux de culture prepares en tubes	19
1.4. Media in tubes for diagnosis of mycobacteria/Milieux de culture en tubes pour diagnostic des mycobactéries	23
2. Dehydrated culture media/milieux de culture deshydrates	24
2.1. Dehydrated culture media/Milieux de culture deshydrates	24
2.2. Ingredients of media/Composants des milieux	39
2.3. Sterile animal blood in bottles/Sang sterile animalier en bouteille	39
2.4. Supplements/Supplements	40
3. Antibiotic discs/disques antibiogramme	43
3.1. Antibiotic paper discs/Disques antibiogramme	43
4. Identification/ d'identification	45
4.1. Diagnostic discs/Disques d'identification	45
4.2. Latex tests/Tests d'agglutination latex	46
4.3. Reagents and tests/Reactifs et tests	46
4.4. Stains/Colorains	47
5. Cryogenic storage of microorganisms - viabank/stockage cryogénique des microorganismes – viabank	47
6. Cell culture media and buffers/milieux de culture de cellules et buffers	47

Prepared culture media

Milieux de culture préparés

Prepared media are manufactured in our production facility, equipped with modern production lines. This provides our company with a high production efficiency and ensures product quality which conforms with the requirements established in the European legislation. Our Quality System is certified to ISO 9001 and ISO 13485 for the design and manufacture of medical devices. The culture media meet the standards of Directive 98/79/CE (1998/10/27) related to medical devices for in vitro diagnostic use. Each lot of prepared media is tested with control strains for growth and (when required) for inhibition. Lot that meets all standards receive The Quality Control Certificate.

Les milieux de cultures préparées sont fabriqués par nos soins à l'aide d'une installation équipée des lignées de production modernes. Cela nous assure une production de haute performance et des produits de qualité conforme aux normes établies par la législation européenne. Notre société est certifiée ISO 9001 et ISO 13485 pour la conception et la fabrication de dispositifs médicaux. Notre système qualité est certifié ISO 9001 et ISO 13485 pour la conception et la fabrication de dispositifs médicaux. Les milieux de culture sont conformes aux normes de la Directive 98/79 / CE relative aux dispositifs médicaux de diagnostic „in vitro”. Chaque lot de milieux de culture ainsi préparé est testé avec des souches de référence de la croissance et (si nécessaire) de l'inhibition.

Prepared culture media on petri dishes

Milieux de culture préparés petri

Description	Cat.No	Package size
Package size 10 or 100 pcs		
Petri dishes Ø 90 mm		
Actinomycete LAB-AGAR™	PP 0122	10 plates
Azide Blood LAB-AGAR™ + 5% SB	PP 1111	10 plates
Azide Blood LAB-AGAR™ + 5% SB + CV	PP 0208	10 plates
BBE LAB-AGAR™	PP 0081	10 plates
BGA LAB-AGAR™ acc. to ISO 6579	PP 1360	10 plates
BHI LAB-AGAR™ + 10 mg/l Teicoplanine	PP 0044	10 plates
BHI LAB-AGAR™ + 2% NaCl	PP 0106	10 plates
BHI LAB-AGAR™ + 2000 mg/l Streptomycin	PP 0009	10 plates
BHI LAB-AGAR™ + 4 mg/l Vancomycin	PP 0010	10 plates
BHI LAB-AGAR™ + 5% SB	PP 0032	10 plates
BHI LAB-AGAR™ + 500 mg/l Gentamicin	PP 0008	10 plates
BHI LAB-AGAR™ + 6 mg/l Vancomycin	PP 0011	10 plates
BHI LAB-AGAR™ + Potassium Tellurite	PP 0125	10 plates
Brain Heart Infusion LAB-AGAR™	PP 0007	10 plates
Bile Chrysoidine Glycerol LAB-AGAR™	PP 0117	10 plates
Bile Esculin Azide LAB-AGAR™	PP 1075	10 plates
Bile Esculin LAB-AGAR™	PP 0019	10 plates
Bismuth Sulfite LAB-AGAR™	PP 1040	10 plates
Blood LAB-AGAR™ + 5% SB	PP 1110	10 plates
Blood LAB-AGAR™ + 7% SB	PP 0195	10 plates
Blood No 2 LAB-AGAR™ + 5% SB	PP 0015	10 plates
Brucella LAB-AGAR™ + 3% SB + wit. K + Hemine	PP 0132	10 plates
Brucella LAB-AGAR™ + 5% HB	PP 0120	10 plates
Brucella LAB-AGAR™ + 5% SB + wit. K + Hemine	PP 0082	10 plates
Brucella LAB-AGAR™ Base	PP 0209	10 plates
Burkholderia cepacia LAB-AGAR™	PP 0066	10 plates
Campylobacter Preston LAB-AGAR™	PP 2160	10 plates
Campylobacter Selective LAB-AGAR™	PP 0048	10 plates
Campylobacter Selective LAB-AGAR™ + 10 % SB	PP 0068	10 plates
Campylobacter Skirrow LAB-AGAR™	PP 0190	10 plates
CCDA LAB-AGAR™	PP 0017	10 plates

Prepared culture media on petri dishes

Milieux de culture préparés petri

Description	Cat.No	Package size
CDC Anaerobe LAB-AGAR™ + 5% SB	PP 0162	10 plates
Cetrimide LAB-AGAR™	PP 1310	10 plates
Charcoal LAB-AGAR™ + 10% HB + Cephalixin	PP 0189	10 plates
Chocolate 2 LAB-AGAR™ + PV	PP 0053	10 plates
Chocolate 3 LAB-AGAR™ + PV	PP 0075	10 plates
Chocolate LAB-AGAR™ + HB + PV	PP 0170	10 plates
Chocolate LAB-AGAR™ + PV	PP 1080	10 plates
Chocolate LAB-AGAR™ + PV + ABV	PP 0210	10 plates
Chocolate LAB-AGAR™ + PV + Bacitracin	PP 1261	10 plates
Chocolate LAB-AGAR™ + PV + VCAT	PP 1082	10 plates
Chocolate LAB-AGAR™ + PV + VCN	PP 1083	10 plates
Chocolate LAB-AGAR™ + PV + VCNT	PP 0200	10 plates
Christensen LAB-GAR™	PP 0214	10 plates
Chromogenic Acinetobacter LAB-AGAR™	PP 2100	10 plates
Chromogenic Acinetobacter MDR LAB-AGAR™	PP 0139	10 plates
Chromogenic C. difficile LAB-AGAR™	PP 0137	10 plates
Chromogenic C. perfringens LAB-AGAR™	PP 0201	10 plates
Chromogenic C3GR LAB-AGAR™	PP 0202	10 plates
Chromogenic Campylobacter LAB-AGAR™	PP 0163	10 plates
Chromogenic Candida LAB-AGAR™	PP 0001B	10 plates
Chromogenic E. coli O157 LAB-AGAR™	PP 0041	10 plates
Chromogenic E. coli STEC LAB-AGAR™	PP 0143	10 plates
Chromogenic E. faecium LAB-AGAR™	PP 0205	10 plates
Chromogenic Enterococcus LAB-AGAR™	PP 0140	10 plates
Chromogenic ESBL LAB-AGAR™	PP 0072	10 plates
Chromogenic ESBL Mod. LAB-AGAR™	PP 0155	10 plates
Chromogenic KPC LAB-AGAR™	PP 0094	10 plates
Chromogenic KPC Mod. LAB-AGAR™	PP 0157	10 plates
Chromogenic Listeria Confirmatory LAB-AGAR™	PP 0156	10 plates
Chromogenic Listeria LAB-AGAR™	PP 0169	10 plates
Chromogenic Malassezia LAB-AGAR™	PP 0203	10 plates
Chromogenic MRSA Modified LAB-AGAR™	PP 0043	10 plates
Chromogenic Orientation LAB-AGAR™	PP 0206	10 plates
Chromogenic Pseudomonas LAB-AGAR™	PP 0097	10 plates
Chromogenic S. aureus Modified LAB-AGAR™	PP 0045	10 plates
Chromogenic Salmonella Plus LAB-AGAR™	PP 0028	10 plates
Chromogenic Salmonella Selective LAB-AGAR™	PP 0115	10 plates
Chromogenic Strepto B LAB-AGAR™	PP 0104	10 plates
Chromogenic Super Carba LAB-AGAR™	PP 0198	10 plates
Chromogenic URI-Colour LAB-AGAR™	PP 0035	10 plates
Chromogenic URI-Colour Modified LAB-AGAR™	PP 0129	10 plates
Chromogenic UTI LAB-AGAR™	PP 0204	10 plates
Chromogenic VRE LAB-AGAR™	PP 0074	10 plates
Chromogenic Yersinia LAB-AGAR™	PP 0144	10 plates
Cled LAB-AGAR™	PP 1030	10 plates
Clostridium difficile LAB-AGAR™	PP 0004	10 plates
Columbia CNA LAB-AGAR™ + 5% SB	PP 1191	10 plates
Columbia LAB-AGAR™ + 10% SB	PP 0065	10 plates
Columbia LAB-AGAR™ + 5 % HB	PP 0020	10 plates
Columbia LAB-AGAR™ + 5% SB	PP 1190	10 plates

Prepared culture media on petri dishes

Milieux de culture préparés petri

Description	Cat.No	Package size
Columbia LAB-AGAR™ + 5% SB + CV	PP 0211	10 plates
Columbia LAB-AGAR™ + 5% SB + V	PP 0212	10 plates
Columbia LAB-AGAR™ + 5% SB + Esculin	PP 0158	10 plates
Columbia LAB-AGAR™ + 5% SB + Potassium Tellurite	PP 0055	10 plates
Columbia LAB-AGAR™ + 7% SB	PP 0059	10 plates
Columbia LAB-AGAR™ + 7% SB + Esculin	PP 0154	10 plates
Columbia LAB-AGAR™ w/o Blood	PP 0213	10 plates
Corynebacterium LAB-AGAR™	PP 0067	10 plates
CR-MOX LAB-AGAR™	PP 0176	10 plates
C-T SMAC LAB-AGAR™	PP 0042	10 plates
Czapek - Dox LAB-AGAR™	PP 0038	10 plates
DCA LAB-AGAR™	PP 0196	10 plates
Dermatophytes LAB-AGAR™	PP 0051	10 plates
DNAse Test LAB-AGAR™	PP 0076	10 plates
Enrichment LAB-AGAR™	PP 1109	10 plates
Enterococcosel LAB-AGAR™ + Potassium Tellurite	PP 0214	10 plates
Enterococcosel LAB-AGAR™ + V	PP 0215	10 plates
Eosin Methylene Blue LAB-AGAR™	PP 1011	10 plates
Esculin Blood LAB-AGAR™	PP 1121	10 plates
Esculin Blood LAB-AGAR™ + 5% SB	PP 0159	10 plates
ESIA LAB-AGAR™ (E. sakazaki)	PP 0031	10 plates
Eugon L.T. LAB-AGAR™	PP 0108	10 plates
Fungisel LAB-AGAR™	PP 0216	10 plates
Gardnerella vaginalis LAB-AGAR™	PP 0005	10 plates
Haemophilus Test Medium	PP 1260	10 plates
Hektoen LAB-AGAR™	PP 1060	10 plates
Helicobacter pylori LAB-AGAR™	PP 0024	10 plates
Karmali LAB-AGAR™	PP 0022	10 plates
Karmali w. Amphotericin LAB-AGAR™	PP 0174	10 plates
Kimming LAB-AGAR™	PP 0147	10 plates
King A LAB-AGARE™ + CN	PP 0070	10 plates
King A LAB-AGAR™	PP 0069	10 plates
King B LAB-AGAR™	PP 0071	10 plates
Legionella BCYE + AB LAB-AGAR™	PP 0178	10 plates
Legionella BCYE + Cephalexin LAB-AGAR™	PP 0188	10 plates
Legionella BCYE LAB-AGAR™	PP 0026	10 plates
Legionella GVPC LAB-AGAR™	PP 0025	10 plates
Legionella MWY LAB-AGAR™	PP 0179	10 plates
Legionella MWY LAB-AGAR™	PP 0179	10 plates
Legionella without Cysteine LAB-AGAR™	PP 0027	10 plates
Mac Conkey LAB-AGAR™	PP 1017	10 plates
Mannitol Salt LAB-AGAR™	PP 1050	10 plates
Mannitol Salt LAB-AGAR™ + Egg Yolk	PP 0217	10 plates
Mannitol Salt LAB-AGAR™ + Oxacillin	PP 0127	10 plates
MRS LAB - AGAR™	PP 0023	10 plates
Mueller Hinton 2 LAB-AGAR™	PP 1170	10 plates
Mueller Hinton 2 LAB-AGAR™ + 2% NaCl	PP 0077	10 plates
Mueller Hinton 2 LAB-AGAR™ + 4% NaCl	PP 0218	10 plates
Mueller Hinton 2 LAB-AGAR™ + 5% SB	PP 1172	10 plates
Mueller Hinton 2 LAB-AGAR™ + Cloxacillin	PP 0073	10 plates

Prepared culture media on petri dishes

Milieux de culture préparés petri

Description	Cat.No	Package size
Mueller Hinton 2 LAB-AGAR™ + NAD + 5% HB	PP 0083	10 plates
Mueller Hinton LAB-AGAR™ + 4% NaCl + Oxacillin	PP 0006	10 plates
Mueller Hinton LAB-AGAR™ + Glucose + Methylene Blue	PP 0105	10 plates
Mycoplasma Selective LAB-AGAR™	PP 0121	10 plates
Nutrient No 1 LAB-AGAR™	PP 0119	10 plates
ORSIM LAB-AGAR™	PP 0128	10 plates
Rice with Tween 80 LAB-AGAR™	PP 0197	10 plates
RPMI MOPS LAB-AGAR™	PP 0095	10 plates
Sabhi LAB-AGAR™	PP 0219	10 plates
Sabouraud 2 % Dextrose with Chloramphenicol LAB-AGAR™	PP 0062	10 plates
Sabouraud 4 % Dextrose with Chloramphenicol LAB-AGAR™	PP 0084	10 plates
Sabouraud Dextrose LAB-AGAR™	PP 1230	10 plates
Sabouraud Dextrose w. Actidione LAB-AGAR™	PP 0220	10 plates
Sabouraud Dextrose with Chloramphenicol and Cycloheximide LAB-AGAR™	PP 0021	10 plates
Sabouraud Dextrose with Chloramphenicol and Gentamicin LAB-AGAR™	PP 1232	10 plates
Sabouraud Dextrose with Chloramphenicol LAB-AGAR™	PP 1231	10 plates
Salmonella Shigella LAB-AGAR™	PP 1250	10 plates
Schaedler Anaerobe LAB-AGAR™	PP 0002	10 plates
Schaedler Anaerobe LAB-AGAR™ + 5% SB	PP 0012	10 plates
Schaedler Anaerobe LAB-AGAR™ + 5% SB + wit. K3	PP 0003	10 plates
Schaedler LAB-AGAR™ + 5% SB + CNA	PP 0039	10 plates
Schaedler LAB-AGAR™ + 5% SB + Gentamicin	PP 0131	10 plates
Schaedler LAB-AGAR™ + 5% SB + K + VA	PP 0014	10 plates
Schaedler LAB-AGAR™ + 5% SB + VA + NEO	PP 0013	10 plates
Sorbitol Mac Conkey LAB-AGAR™	PP 1021	10 plates
Streptococcus Selective LAB-AGAR™	PP 0146	10 plates
TCBS LAB-AGAR™	PP 0036	10 plates
Trypticasein Soy LAB-AGAR™	PP 1180	10 plates
Trypticasein Soy LAB-AGAR™ + 5% SB	PP 1181	10 plates
TYCSB Agar	PP 0225	10 plates
Wilkins Chalgren LAB-AGAR™	PP 0221	10 plates
Wilkins Chalgren LAB-AGAR™ + 5% HB	PP 0222	10 plates
Wilkins Chalgren LAB-AGAR™ + 5% SB	PP 0063	10 plates
XLD LAB-AGAR™	PP 1330	10 plates
Yeast Nitrogen Modified LAB-AGAR™	PP 0052	10 plates
Yersinia CIN LAB-AGAR™	PP 1090	10 plates
Petri dishes Ø 140 mm		
Mueller Hinton 2 LAB-AGAR™ (140 mm)	PP 0016	10 plates
Mueller Hinton 2 LAB-AGAR™ + 5% SB (140 mm)	PP 0018	10 plates
Mueller Hinton 2 LAB-AGAR™ + NAD + 5% HB (140 mm)	PP 0092	10 plates
RPMI MOPS LAB-AGAR™ (140 mm)	PP 0124	10 plates
Two-sections Petri dishes Ø 90 mm		
BHI LAB-AGAR™ + 500 mg/l Gentamicin / BHI LAB-AGAR™ + 2000 mg/l Streptomycin	PD 007	10 plates
BHI LAB-AGAR™ + Potassium Tellurite / BHI LAB-AGAR™ + Potassium Tellurite	PD 052	10 plates
Bismuth Sulfite LAB-AGAR™ / Bismuth Sulfite LAB-AGAR™	PD 055	10 plates
Blood LAB-AGAR™ + 5% SB / MacConkey LAB-AGAR™	PD 001	10 plates
Chocolate LAB-AGAR™ + PV / Chocolate LAB-AGAR™ + PV + Bacitracin	PD 049	10 plates
Chromogenic Candida LAB-AGAR™ / Chromogenic Candida LAB-AGAR™	PD 026	10 plates
Chromogenic Candida LAB-AGAR™ / Sabouraud Dext. LAB-AGAR + Ge + C	PD 017	10 plates
Chromogenic ESBL LAB-AGAR™ / Chromogenic KPC LAB-AGAR™	PD 062	10 plates

Prepared culture media on petri dishes

Milieux de culture préparés petri

Description	Cat.No	Package size
Chromogenic ESBL LAB-AGAR™ / Chromogenic Orientation LAB-AGAR™	PD 063	10 plates
Chromogenic ESBL LAB-AGAR™ / Chromogenic VRE LAB-AGAR™	PD 042	10 plates
Chromogenic Strepto B / Chromogenic Listeria LAB-AGAR™	PD 045	10 plates
Chromogenic S. aureus LAB-AGAR™ / Chromogenic MRSA LAB-AGAR™	PD 041	10 plates
Chromogenic MRSA LAB-AGAR™ / Chromogenic VRE LAB-AGAR™	PD 064	10 plates
Chromogenic Salmonella LAB-AGAR™ / Chromogenic E. coli O157	PD 065	10 plates
Chromogenic Salmonella LAB-AGAR™ / Chromogenic Salmonella LAB-AGAR™	PD 044	10 plates
Chromogenic Salmonella LAB-AGAR™ / Hektoen LAB-AGAR™	PD 024	10 plates
Chromogenic Strepto B LAB-AGAR™ / Columbia CNA LAB-AGAR™ + 5% SB	PD 066	10 plates
Chromogenic Strepto B LAB-AGAR™ / Columbia LAB-AGAR™ + 5% SB	PD 016	10 plates
Chromogenic KPC LAB-AGAR™ / Chromogenic VRE LAB-AGAR™	PD 046	10 plates
Chromogenic Salmonella / XLD LAB-AGAR™	PD 023	10 plates
Chromogenic Strepto B / Chromogenic Strepto B LAB-AGAR™	PD 046	10 plates
Chromogenic Strepto B / Columbia LAB-AGAR™ + 5% SB	PD 016	10 plates
Chromogenic URI LAB-AGAR™ / Columbia LAB-AGAR™ + 5% SB	PD 018	10 plates
Chromogenic URI LAB-AGAR™ / Columbia CNA LAB-AGAR™ + 5% SB	PD 025	10 plates
Columbia LAB-AGAR™ / Chocolate LAB-AGAR™ + PV	PD 038	10 plates
Columbia LAB-AGAR™ + 5% SB / Chocolate LAB-AGAR™ + PV + Bacitracin	PD 067	10 plates
Columbia LAB-AGAR™ + 5% SB / Chromogenic Orientation LAB-AGAR™	PD 068	10 plates
Columbia LAB-AGAR™ + 5% SB / Columbia LAB-AGAR™ + 5% SB	PD 032	10 plates
Columbia LAB-AGAR™ + 5% SB / MacConkey LAB-AGAR™	PD 014	10 plates
Columbia LAB-AGAR™ + 5% SB / Mannitol Salt LAB-AGAR™	PD 048	10 plates
Columbia LAB-AGAR™ + 5% SB / Sabouraud Dextrose with Chloramphenicol LAB-AGAR™	PD 031	10 plates
Columbia CNA LAB-AGAR™ + 5% SB + Esculin / Chromogenic URI LAB-AGAR™	PD 030	10 plates
Columbia CNA LAB-AGAR™ + 5% SB / Bile Esculin Azide LAB-AGAR™	PD 008	10 plates
Columbia CNA LAB-AGAR™ + 5% SB / Chocolate LAB-AGAR™ + PV + Bacitracin	PD 010	10 plates
Columbia CNA LAB-AGAR™ + 5% SB / Chocolate LAB-AGAR™ + PV	PD 009	10 plates
Columbia CNA LAB-AGAR™ + 5% SB / Chromogenic Candida LAB-AGAR™	PD 069	10 plates
Columbia CNA LAB-AGAR™ + 5% KB / Chromogenic Orientation LAB-AGAR™	PD 070	10 plates
Columbia CNA LAB-AGAR™ + 5% KB / MacConkey LAB-AGAR™	PD 011	10 plates
Columbia CNA LAB-AGAR™ + 5% KB / Mannitol Salt LAB-AGAR™	PD 012	10 plates
Columbia CNA LAB-AGAR™ + 5% SB / Sabouraud with Chloramphenicol and Gentamicin LAB-AGAR™	PD 013	10 plates
Dermatophytes LAB-AGAR™ / Dermatophytes LAB-AGAR™	PD 035	10 plates
Enterococcosel LAB-AGAR™ / Enterococcosel LAB-AGAR™ + Potassium Tellurite	PD 071	10 plates
Eosin Methylene Blue LAB-AGAR™ / Blood LAB-AGAR™ + 5% SB	PD 002	10 plates
Eosin Methylene Blue LAB-AGAR™ / Mannitol Salt LAB-AGAR™	PD 003	10 plates
Legionella BCYE LAB-AGAR™ / Legionella w/o Cysteine LAB-AGAR™	PD 015	10 plates
Legionella BCYE LAB-AGAR™ / Legionella BCYE LAB-AGAR™	PD 072	10 plates
Mac Conkey Sorbitol LAB-AGAR™ / Yersinia CIN LAB-AGAR™	PD 073	10 plates
Mac Conkey LAB-AGAR™ / Chromogenic Candida LAB-AGAR™	PD 061	10 plates
Mac Conkey LAB-AGAR™ / CLED LAB-AGAR™	PD 060	10 plates
Mac Conkey LAB-AGAR™ / Mac Conkey LAB-AGAR™	PD 053	10 plates
MacConkey LAB-AGAR™ / Sabouraud Dextrose with Chloramphenicol LAB-AGAR™	PD 029	10 plates
Mac Conkey LAB-AGAR™ / Sabouraud with Chloramphenicol and Gentamicin LAB-AGAR™	PD 039	10 plates
Mueller Hinton 2 LAB-AGAR™ + 5% SB / Mueller Hinton 2 LAB-AGAR™ + 5% SB	PD 074	10 plates
Mueller Hinton 2 LAB-AGAR™ / Mueller Hinton 2 LAB-AGAR™	PD 033	10 plates
ORSIM LAB-AGAR™ / ORSIM LAB-AGAR™	PD 027	10 plates
Sabouraud LAB-AGAR™ + Chloramphenicol + Actidione / Fungisel LAB-AGAR™	PD 075	10 plates
Sabouraud Dextrose with Chloramphenicol LAB-AGAR™ / Dermatophytes LAB-AGAR™	PD 005	10 plates

Prepared culture media on petri dishes

Milieux de culture préparés petri

Description	Cat.No	Package size
Sabouraud LAB-AGAR™ + Chloramphenicol / Fungisel LAB-AGAR™	PD 076	10 plates
Sabouraud Dext. w. Chloramph. LAB-AGAR™ / Sabouraud Dext. w. Actidione Lab-AGAR™	PD 043	10 plates
Sabouraud LAB-AGAR™ / Chromogenic Candida LAB-AGAR™	PD 077	10 plates
Sabouraud Dextrose LAB-AGAR™ / Dermatophytes LAB-AGAR™	PD 004	10 plates
Sabouraud LAB-AGAR™ / Fungisel LAB-AGAR™	PD 050	10 plates
Sabouraud Dext. LAB-AGAR™ / Sabouraud Dex. w. Actidione Lab-AGAR™	PD 036	10 plates
Sabouraud w. Chloramphenicol + Actidione LAB-AGAR™ / Sabouraud w. Chloramphenicol + Actidione + Phenol Red LAB-AGAR™	PD 057	10 plates
Salmonella Shigella LAB-AGAR™ / Salmonella Shigella LAB-AGAR™	PD 054	10 plates
Schaedler LAB-AGAR™ + 5% SB / Schaedler LAB-AGAR™ + 5% SB + NEO + VA	PD 022	10 plates
Schaedler LAB-AGAR™ + 5% SB / Schaedler LAB-AGAR™ + 5% SB + K + VA	PD 021	10 plates
Sorbitol Mac Conkey LAB-AGAR™ / Sorbitol Mac Conkey LAB-AGAR™	PD 056	10 plates
SS LAB-AGAR™ / Bismuth Sulfite LAB-AGAR™	PD 006	10 plates
SS LAB-AGAR™ / Chromogenic Salmonella LAB-AGAR™	PD 020	10 plates
SS LAB-AGAR™ / Hektoen LAB-AGAR™	PD 019	10 plates
SS LAB-AGAR™ / XLD LAB-AGAR™	PD 037	10 plates
Wilkins Chalgren LAB-AGAR™ + 5% SB / Wilkins Chalgren LAB-AGAR™ + 5% SB	PD 034	10 plates
Three-sections Petri dishes Ø 90 mm		
available on request		
Four-sections Petri dishes Ø 90 mm		
available on request		

Control of surfaces Petri dishes Ø 55 mm

Total bacteria count

Description	Cat.No	Package size
TSA IRR. LAB-AGAR™	BRD 10	4 x 5 pcs
TSA IRR. LAB-AGAR™ + Letheen	BRD 11	4 x 5 pcs
TSA IRR. LAB-AGAR™ + Letheen + Tween	BRD 12	4 x 5 pcs
TSA IRR. LAB-AGAR™ + 3 N	BRD 19	4 x 5 pcs
TSA IRR. LAB-AGAR™ + 4 N	BRD 14	4 x 5 pcs
Yeast and mold irradiated		
Sabouraud Dextrose IRR. LAB-AGAR™	BRD 04	4 x 5 pcs
Sabouraud Dextrose IRR. LAB-AGAR™ + Letheen	BRD 05	4 x 5 pcs
Sabouraud Dextrose IRR. LAB-AGAR™ + Letheen + Tween	BRD 06	4 x 5 pcs
Sabouraud Dextrose IRR. LAB-AGAR™ + 4N	BRD 26	4 x 5 pcs
Sabouraud Dextrose IRR. LAB-AGAR™ + 4N	BRD 26B	4 x 5 pcs
Sabouraud Dextrose with Chloramphenicol IRR. LAB-AGAR™	BRD 07	4 x 5 pcs
Sabouraud Dextrose with Chloramphenicol IRR. LAB-AGAR™ + Letheen	BRD 08	4 x 5 pcs
Sabouraud Dextrose with Chloramphenicol IRR. LAB-AGAR™ + Letheen + Tween	BRD 09	4 x 5 pcs
Enterobacteriaceae		
VRBG LAB-AGAR™ + LETHEEN + TWEEN	BR 13	4 x 5 pcs
VRBG LAB-AGAR™	BR 17	4 x 5 pcs
VRBL LAB-AGAR™	BR 18	4 x 5 pcs
VRBL LAB-AGAR™ + 4N	BR 32	4 x 5 pcs
PSEUDOMONAS SPP		
Cetrimide LAB-AGAR™	BR 16	4 x 5 pcs

Control of surfaces Petri dishes Ø 55 mm

Total bacteria count

Description	Cat.No	Package size
Chromogenic Pseudomonas LAB-AGAR™	BR 51	4 x 5 pcs
Cetrimide LAB-AGAR™ + Letheen + Tween	BR 20	4 x 5 pcs
Staphylococcus spp.		
BAIRD PARKER LAB-AGAR™ + LETHEEN + TWEEN	BR 28	4 x 5 pcs
MANNITOL SALT LAB- AGAR™	BR 15	4 x 5 pcs
MANNITOL SALT LAB- AGAR™ + LETHEEN + TWEEN	BR 21	4 x 5 pcs
Salmonella spp.		
Chromogenic Salmonella Selective LAB-AGAR™	BR 25	4 x 5 pcs
Chromogenic Salmonella Selective LAB-AGAR™ + Letheen + Tween	BR 24	4 x 5 pcs
Listeria spp.		
Chromogenic Listeria LAB-AGAR™	BR 31	4 x 5 pcs

Prepared culture media in bottles

Milieux de culture prepares en bouteilles

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
0,85% Sodium Chloride	BT 5054.01	6 x 100 ml
0,85% Sodium Chloride	BT 5054.02	6 x 200 ml
0,85% Sodium Chloride	BT 5054.05	6 x 500 ml
0,9% Sodium Chloride	BT 5045.01	6 x 100 ml
0,9% Sodium Chloride	BT 5045.02	6 x 200 ml
0,9% Sodium Chloride	BT 5045.05	6 x 500 ml
1% Glucose Broth	BT 5064.01	6 x 100 ml
1% Glucose Broth	BT 5064.02	6 x 200 ml
1% Glucose Broth	BT 5064.05	6 x 500 ml
10% Lactose Broth	BT 5052.01	6 x 100 ml
10% Lactose Broth	BT 5052.02	6 x 200 ml
10% Lactose Broth	BT 5052.05	6 x 500 ml
Acetamide Medium	BT 5214.01	6 x 100 ml
Acetamide Medium	BT 5214.02	6 x 200 ml
Acetamide Medium	BT 5214.05	6 x 500 ml
Actinomycete LAB-AGAR™	BT 5181.01	6 x 100 ml
Actinomycete LAB-AGAR™	BT 5181.02	6 x 200 ml
Actinomycete LAB-AGAR™	BT 5181.05	6 x 500 ml
Alkaline Water	BT 5262	6 x 225 ml
Alkaline Water 10x	BT 5041.01	6 x 100 ml
Alkaline Water 10x	BT 5041.02	6 x 200 ml
Alkaline Water 10x	BT 5041.05	6 x 500 ml
Alkaline Peptone Water	BT 5273.01	6 x 100 ml
Alkaline Peptone Water	BT 5273.02	6 x 200 ml
Alkaline Peptone Water	BT 5273.05	6 x 500 ml
Andrade Medium Base	BT 5235.01	6 x 100 ml
Andrade Medium Base	BT 5235.02	6 x 200 ml
Andrade Medium Base	BT 5235.05	6 x 500 ml
Arabinose Broth	BT 5053.01	6 x 100 ml
Arabinose Broth	BT 5053.02	6 x 200 ml
Arabinose Broth	BT 5053.05	6 x 500 ml
Arginine acc. to Falkow Broth	BT 5049.01	6 x 100 ml
Arginine acc. to Falkow Broth	BT 5049.02	6 x 200 ml
Arginine acc. to Falkow Broth	BT 5049.05	6 x 500 ml
Azide Dextrose Broth	BT 5270.01	6 x 100 ml
Azide Dextrose Broth	BT 5270.02	6 x 200 ml
Azide Dextrose Broth	BT 5270.05	6 x 500 ml
Azide LAB-AGAR™ + CV	BT 5103.01	6 x 100 ml
Azide LAB-AGAR™ + CV	BT 5103.02	6 x 200 ml
Azide LAB-AGAR™ + CV	BT 5103.05	6 x 500 ml
BGA LAB-AGAR™ acc. to ISO 6579	BT 5096.01	6 x 100 ml
BGA LAB-AGAR™ acc. to ISO 6579	BT 5096.02	6 x 200 ml
BGA LAB-AGAR™ acc. to ISO 6579	BT 5096.05	6 x 500 ml
Bile Esculin Azide LAB-AGAR™	BT 5088.01	6 x 100 ml
Bile Esculin Azide LAB-AGAR™	BT 5088.02	6 x 200 ml
Bile Esculin Azide LAB-AGAR™	BT 5088.05	6 x 500 ml
Bile Salts	BT 5271.01	6 x 100 ml

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
Bile Salts	BT 5271.02	6 x 200 ml
Bile Salts	BT 5271.05	6 x 500 ml
Bismuth sulfite LAB-AGAR™	BT 5003.01	6 x 100 ml
Bismuth sulfite LAB-AGAR™	BT 5003.02	6 x 200 ml
Bismuth sulfite LAB-AGAR™	BT 5003.05	6 x 500 ml
Blood LAB-AGAR™ Base	BT 5034.01	6 x 100 ml
Blood LAB-AGAR™ Base	BT 5034.02	6 x 200 ml
Blood LAB-AGAR™ Base	BT 5034.05	6 x 500 ml
Brain Heart Infusion Broth	BT 5008.01	6 x 100 ml
Brain Heart Infusion Broth	BT 5008.02	6 x 200 ml
Brain Heart Infusion Broth	BT 5008.05	6 x 500 ml
Brain Heart Infusion Broth	BT 5197	6 x 30 ml
Brain Heart Infusion Broth	BT 5226	6 x 70 ml
Brain Heart Infusion LAB-AGAR™	BT 5009.01	6 x 100 ml
Brain Heart Infusion LAB-AGAR™	BT 5009.02	6 x 200 ml
Brain Heart Infusion LAB-AGAR™	BT 5009.05	6 x 500 ml
Branhamella LAB-AGAR™	BT 5272.01	6 x 100 ml
Branhamella LAB-AGAR™	BT 5272.02	6 x 200 ml
Branhamella LAB-AGAR™	BT 5272.05	6 x 500 ml
Brucella LAB-AGAR™ Base	BT 5274.01	6 x 100 ml
Brucella LAB-AGAR™ Base	BT 5274.02	6 x 200 ml
Brucella LAB-AGAR™ Base	BT 5274.05	6 x 500 ml
Buffered Salt Dilution (PBS)	BT 5105.01	6 x 100 ml
Buffered Salt Dilution (PBS)	BT 5105.02	6 x 200 ml
Buffered Salt Dilution (PBS)	BT 5105.05	6 x 500 ml
Bulion (Broth)+ 6,5% NaCl	BT 5079.01	6 x 100 ml
Bulion (Broth)+ 6,5% NaCl	BT 5079.02	6 x 200 ml
Bulion (Broth)+ 6,5% NaCl	BT 5079.05	6 x 500 ml
Bulion (Broth) pH 9,6	BT 5081.01	6 x 100 ml
Bulion (Broth) pH 9,6	BT 5081.02	6 x 200 ml
Bulion (Broth) pH 9,6	BT 5081.05	6 x 500 ml
Bulion zwykły (Nutritious broth)	BT 5084.01	6 x 100 ml
Bulion zwykły (Nutritious broth)	BT 5084.02	6 x 200 ml
Bulion zwykły (Nutritious broth)	BT 5084.05	6 x 500 ml
Cary-Blair Transport Medium	BT 5237.01	6 x 100 ml
Cary-Blair Transport Medium	BT 5237.02	6 x 200 ml
Cary-Blair Transport Medium	BT 5237.05	6 x 500 ml
Celobiose Broth	BT 5085.01	6 x 100 ml
Celobiose Broth	BT 5085.02	6 x 200 ml
Celobiose Broth	BT 5085.05	6 x 500 ml
Cetrimide LAB-AGAR™	BT 5122.01	6 x 100 ml
Cetrimide LAB-AGAR™	BT 5122.02	6 x 200 ml
Cetrimide LAB-AGAR™	BT 5122.05	6 x 500 ml
Cetrimide LAB-AGAR™	BT 7004	6 x 250 ml
Chocolate LAB-AGAR™ Base	BT 5017.01	6 x 100 ml
Chocolate LAB-AGAR™ Base	BT 5017.02	6 x 200 ml
Chocolate LAB-AGAR™ Base	BT 5017.05	6 x 500 ml
Chromogenic Candida LAB-AGAR™	BT 5222.01	6 x 100 ml
Chromogenic Candida LAB-AGAR™	BT 5222.02	6 x 200 ml
Chromogenic Candida LAB-AGAR™	BT 5222.05	6 x 500 ml

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
Chromogenic Candida LAB-AGAR™	BT 7001	6 x 250 ml
Chromogenic E. coli O157 LAB-AGAR™	BT 5071.01	6 x 100 ml
Chromogenic E. coli O157 LAB-AGAR™	BT 5071.02	6 x 200 ml
Chromogenic E. coli O157 LAB-AGAR™	BT 5071.05	6 x 500 ml
Chromogenic Orientation LAB-AGAR™	BT 5276.01	6 x 100 ml
Chromogenic Orientation LAB-AGAR™	BT 5276.02	6 x 200 ml
Chromogenic Orientation LAB-AGAR™	BT 5276.05	6 x 500 ml
Chromogenic Pseudomonas LAB-AGAR™	BT 5277.01	6 x 100 ml
Chromogenic Pseudomonas LAB-AGAR™	BT 5277.02	6 x 200 ml
Chromogenic Pseudomonas LAB-AGAR™	BT 5277.05	6 x 500 ml
Chromogenic Salmonella Plus LAB-AGAR™	BT 5275.01	6 x 100 ml
Chromogenic Salmonella Plus LAB-AGAR™	BT 5275.02	6 x 200 ml
Chromogenic Salmonella Plus LAB-AGAR™	BT 5275.05	6 x 500 ml
Chromogenic Salmonella Selective LAB-AGAR™	BT 5072.01	6 x 100 ml
Chromogenic Salmonella Selective LAB-AGAR™	BT 5072.02	6 x 200 ml
Chromogenic Salmonella Selective LAB-AGAR™	BT 5072.05	6 x 500 ml
Chromogenic Staph. Aureus LAB-AGAR™	BT 5278.01	6 x 100 ml
Chromogenic Staph. Aureus LAB-AGAR™	BT 5278.02	6 x 200 ml
Chromogenic Staph. Aureus LAB-AGAR™	BT 5278.05	6 x 500 ml
Chromogenic URI-COLOUR LAB-AGAR™	BT 5073.01	6 x 100 ml
Chromogenic URI-COLOUR LAB-AGAR™	BT 5073.02	6 x 200 ml
Chromogenic URI-COLOUR LAB-AGAR™	BT 5073.05	6 x 500 ml
Citrate Christensen LAB-AGAR™	BT 5042.01	6 x 100 ml
Citrate Christensen LAB-AGAR™	BT 5042.02	6 x 200 ml
Citrate Christensen LAB-AGAR™	BT 5042.05	6 x 500 ml
Christensen LAB-AGAR™	BT 5279.01	6 x 100 ml
Christensen LAB-AGAR™	BT 5279.02	6 x 200 ml
Christensen LAB-AGAR™	BT 5279.05	6 x 500 ml
Cled LAB-AGAR™	BT 5102.01	6 x 100 ml
Cled LAB-AGAR™	BT 5102.02	6 x 200 ml
Cled LAB-AGAR™	BT 5102.05	6 x 500 ml
Columbia LAB-AGAR™ Base	BT 5011.01	6 x 100 ml
Columbia LAB-AGAR™ Base	BT 5011.02	6 x 200 ml
Columbia LAB-AGAR™ Base	BT 5011.05	6 x 500 ml
Columbia CNA LAB-AGAR™ Base	BT 5280.01	6 x 100 ml
Columbia CNA LAB-AGAR™ Base	BT 5280.02	6 x 200 ml
Columbia CNA LAB-AGAR™ Base	BT 5280.05	6 x 500 ml
DCA LAB-AGAR™	BT 5281.01	6 x 100 ml
DCA LAB-AGAR™	BT 5281.02	6 x 200 ml
DCA LAB-AGAR™	BT 5281.05	6 x 500 ml
Distilled Water	BT 5172.01	6 x 100 ml
Distilled Water	BT 5172.02	6 x 200 ml
Distilled Water	BT 5172.05	6 x 500 ml
DNase Test LAB-AGAR™	BT 5078.01	6 x 100 ml
DNase Test LAB-AGAR™	BT 5078.02	6 x 200 ml
DNase Test LAB-AGAR™	BT 5078.05	6 x 500 ml
Drigalski LAB-AGAR™	BT 5282.01	6 x 100 ml
Drigalski LAB-AGAR™	BT 5282.02	6 x 200 ml
Drigalski LAB-AGAR™	BT 5282.05	6 x 500 ml
Dulcitol Broth	BT 5043.01	6 x 100 ml

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
Dulcitol Broth	BT 5043.02	6 x 200 ml
Dulcitol Broth	BT 5043.05	6 x 500 ml
Dermatophytes LAB-AGAR™	BT 5282.01	6 x 100 ml
Dermatophytes LAB-AGAR™	BT 5282.02	6 x 200 ml
Dermatophytes LAB-AGAR™	BT 5282.05	6 x 500 ml
Enrichment Broth	BT 5283.01	6 x 100 ml
Enrichment Broth	BT 5283.02	6 x 200 ml
Enrichment Broth	BT 5283.05	6 x 500 ml
Enrichment LAB-AGAR™	BT 5284.01	6 x 100 ml
Enrichment LAB-AGAR™	BT 5284.02	6 x 200 ml
Enrichment LAB-AGAR™	BT 5284.05	6 x 500 ml
Enterococcus Arabinose Broth	BT 5062.01	6 x 100 ml
Enterococcus Arabinose Broth	BT 5062.02	6 x 200 ml
Enterococcus Arabinose Broth	BT 5062.05	6 x 500 ml
Eosin Methylene Blue LAB-AGAR™	BT 5001.01	6 x 100 ml
Eosin Methylene Blue LAB-AGAR™	BT 5001.02	6 x 200 ml
Eosin Methylene Blue LAB-AGAR™	BT 5001.05	6 x 500 ml
Esterase Tween LAB-AGAR™	BT 5230.01	6 x 100 ml
Eugon Broth	BT 5285.01	6 x 100 ml
Eugon Broth	BT 5285.02	6 x 200 ml
Eugon Broth	BT 5285.05	6 x 500 ml
Ewing Malonate Broth Modified	BT 5018.01	6 x 100 ml
Ewing Malonate Broth Modified	BT 5018.02	6 x 200 ml
Ewing Malonate Broth Modified	BT 5018.05	6 x 500 ml
Falkow Broth Base	BT 5051.01	6 x 100 ml
Falkow Broth Base	BT 5051.02	6 x 200 ml
Falkow Broth Base	BT 5051.05	6 x 500 ml
Fungisel LAB-AGAR™	BT 5286.01	6 x 100 ml
Fungisel LAB-AGAR™	BT 5286.02	6 x 200 ml
Fungisel LAB-AGAR™	BT 5286.05	6 x 500 ml
Garda LAB-AGAR™	BT 5076.01	6 x 100 ml
Garda LAB-AGAR™	BT 5076.02	6 x 200 ml
Garda LAB-AGAR™	BT 5076.05	6 x 500 ml
Glycerol Broth	BT 5086.01	6 x 100 ml
Glycerol Broth	BT 5086.02	6 x 200 ml
Glycerol Broth	BT 5086.05	6 x 500 ml
Hektoen LAB-AGAR™	BT 5111.01	6 x 100 ml
Hektoen LAB-AGAR™	BT 5111.02	6 x 200 ml
Hektoen LAB-AGAR™	BT 5111.05	6 x 500 ml
Hugh Leifson Medium	BT 5287.01	6 x 100 ml
Hugh Leifson Medium	BT 5287.02	6 x 200 ml
Hugh Leifson Medium	BT 5287.05	6 x 500 ml
Inositol Broth	BT 5087.01	6 x 100 ml
Inositol Broth	BT 5087.02	6 x 200 ml
Inositol Broth	BT 5087.05	6 x 500 ml
King A LAB-AGAR™	BT 5288.01	6 x 100 ml
King A LAB-AGAR™	BT 5288.02	6 x 200 ml
King A LAB-AGAR™	BT 5288.05	6 x 500 ml
King B LAB-AGAR™	BT 5110.01	6 x 100 ml
King B LAB-AGAR™	BT 5110.02	6 x 200 ml

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
King B LAB-AGAR™	BT 5110.05	6 x 500 ml
Kligler Iron LAB-AGAR™	BT 5006.01	6 x 100 ml
Kligler Iron LAB-AGAR™	BT 5006.02	6 x 200 ml
Kligler Iron LAB-AGAR™	BT 5006.05	6 x 500 ml
Lysine acc. to Falkow Broth	BT 5050.01	6 x 100 ml
Lysine acc. to Falkow Broth	BT 5050.02	6 x 200 ml
Lysine acc. to Falkow Broth	BT 5050.05	6 x 500 ml
Mac Conkey LAB-AGAR™	BT 5010.01	6 x 100 ml
Mac Conkey LAB-AGAR™	BT 5010.02	6 x 200 ml
Mac Conkey LAB-AGAR™	BT 5010.05	6 x 500 ml
Maltose Broth	BT 5080.01	6 x 100 ml
Maltose Broth	BT 5080.02	6 x 200 ml
Maltose Broth	BT 5080.05	6 x 500 ml
Mannitol Broth	BT 5048.01	6 x 100 ml
Mannitol Broth	BT 5048.02	6 x 200 ml
Mannitol Broth	BT 5048.05	6 x 500 ml
Mannitol Salt LAB-AGAR™	BT 5104.01	6 x 100 ml
Mannitol Salt LAB-AGAR™	BT 5104.02	6 x 200 ml
Mannitol Salt LAB-AGAR™	BT 5104.05	6 x 500 ml
Melibiose Broth (ISO10273)(Yersinia)	BT 5243.01	6 x 100 ml
Melibiose Broth (ISO10273)(Yersinia)	BT 5243.02	6 x 200 ml
Melibiose Broth (ISO10273)(Yersinia)	BT 5243.05	6 x 500 ml
Melobiose Broth	BT 5093.01	6 x 100 ml
Melobiose Broth	BT 5093.02	6 x 200 ml
Melobiose Broth	BT 5093.05	6 x 500 ml
Motility LAB-AGAR™	BT 5026.01	6 x 100 ml
Motility LAB-AGAR™	BT 5026.02	6 x 200 ml
Motility LAB-AGAR™	BT 5026.05	6 x 500 ml
MRS LAB - AGAR™	BT 5031.01	6 x 100 ml
MRS LAB - AGAR™	BT 5031.02	6 x 200 ml
MRS LAB - AGAR™	BT 5031.05	6 x 500 ml
MRS Broth	BT 5164.01	6 x 100 ml
MRS Broth	BT 5164.02	6 x 200 ml
MRS Broth	BT 5164.05	6 x 500 ml
MRVP Medium	BT 5036.01	6 x 100 ml
MRVP Medium	BT 5036.02	6 x 200 ml
MRVP Medium	BT 5036.05	6 x 500 ml
Mueller Hinton 2 LAB-AGAR™	BT 5002.01	6 x 100 ml
Mueller Hinton 2 LAB-AGAR™	BT 5002.02	6 x 200 ml
Mueller Hinton 2 LAB-AGAR™	BT 5002.05	6 x 500 ml
Mueller Hinton Broth	BT 5055.01	6 x 100 ml
Mueller Hinton Broth	BT 5055.02	6 x 200 ml
Mueller Hinton Broth	BT 5055.05	6 x 500 ml
Nitrate Broth	BT 5056.01	6 x 100 ml
Nitrate Broth	BT 5056.02	6 x 200 ml
Nitrate Broth	BT 5056.05	6 x 500 ml
Nutrient Broth	BT 7012.01	6 x 100 ml
Nutrient Broth	BT 7012.02	6 x 200 ml
Nutrient Broth	BT 7012.05	6 x 500 ml
Nutrient G Broth	BT 5291.01	6 x 100 ml

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
Nutrient G Broth	BT 5291.02	6 x 200 ml
Nutrient G Broth	BT 5292.05	6 x 500 ml
Nutrient Pneumococcus Broth	BT 5268.01	6 x 100 ml
Nutrient Pneumococcus Broth	BT 5268.02	6 x 200 ml
Nutrient Pneumococcus Broth	BT 5268.05	6 x 500 ml
Nutrient LAB-AGAR™	BT 5126.01	6 x 100 ml
Nutrient LAB-AGAR™	BT 5126.02	6 x 200 ml
Nutrient LAB-AGAR™	BT 5126.05	6 x 500 ml
Nutrient LAB-AGAR™ at pH 7,0	BT 5289.01	6 x 100 ml
Nutrient LAB-AGAR™ at pH 7,0	BT 5289.02	6 x 200 ml
Nutrient LAB-AGAR™ at pH 7,0	BT 5289.05	6 x 500 ml
Nutrient No 1 LAB-AGAR™	BT 5195.01	6 x 100 ml
Nutrient No 1 LAB-AGAR™	BT 5195.02	6 x 200 ml
Nutrient No 1 LAB-AGAR™	BT 5195.05	6 x 500 ml
Nutrient 2 LAB-AGAR™	BT 5075.01	6 x 100 ml
Nutrient 2 LAB-AGAR™	BT 5075.02	6 x 200 ml
Nutrient 2 LAB-AGAR™	BT 5075.05	6 x 500 ml
Nutrient S LAB-AGAR™ (for bacterial strains vitality preservation)	BT 5116.01	6 x 100 ml
Nutrient S LAB-AGAR™ (for bacterial strains vitality preservation)	BT 5116.02	6 x 200 ml
Nutrient S LAB-AGAR™ (for bacterial strains vitality preservation)	BT 5116.05	6 x 500 ml
Ornithine acc. to Falkow Broth	BT 5059.01	6 x 100 ml
Ornithine acc. to Falkow Broth	BT 5059.02	6 x 200 ml
Ornithine acc. to Falkow Broth	BT 5059.05	6 x 500 ml
PBS (Buffered Salt Dilution)	BT 5105.01	6 x 100 ml
PBS (Buffered Salt Dilution)	BT 5105.02	6 x 200 ml
PBS (Buffered Salt Dilution)	BT 5105.05	6 x 500 ml
PBS with calcium and magnesium ions	BT 5290.01	6 x 100 ml
PBS with calcium and magnesium ions	BT 5290.02	6 x 200 ml
PBS with calcium and magnesium ions	BT 5290.05	6 x 500 ml
Peptone Water	BT 5046.01	6 x 100 ml
Peptone Water	BT 5046.01	6 x 200 ml
Peptone Water	BT 5046.01	6 x 500 ml
Phenylalanine LAB-AGAR™	BT 5029.01	6 x 100 ml
Phenylalanine LAB-AGAR™	BT 5029.02	6 x 200 ml
Phenylalanine LAB-AGAR™	BT 5029.05	6 x 500 ml
Phenol Red Broth Base	BT 5229.01	6 x 100 ml
Pyrazinamidaze LAB-AGAR™	BT 5231.01	6 x 100 ml
Potato Dextrose LAB-AGAR™	BT 5136.01	6 x 100 ml
Potato Dextrose LAB-AGAR™	BT 5136.01	6 x 200 ml
Potato Dextrose LAB-AGAR™	BT 5136.01	6 x 500 ml
Preston Campylobacter LAB-AGAR™ Base	BT 5293.01	6 x 100 ml
Preston Campylobacter LAB-AGAR™ Base	BT 5293.02	6 x 200 ml
Preston Campylobacter LAB-AGAR™ Base	BT 5293.05	6 x 500 ml
Rafinose Broth	BT 5090.01	6 x 100 ml
Rafinose Broth	BT 5090.02	6 x 200 ml
Rafinose Broth	BT 5090.05	6 x 500 ml
Rice with Tween 80 LAB-AGAR™	BT 5294.01	6 x 100 ml
Rice with Tween 80 LAB-AGAR™	BT 5294.02	6 x 200 ml
Rice with Tween 80 LAB-AGAR™	BT 5294.05	6 x 500 ml
Ringer Solution	BT 5132	6 x 90 ml

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
Ringer Solution	BT 6000	6 x 180 ml
Ringer Solution	BT 5140.01	6 x 100 ml
Ringer Solution	BT 5140.02	6 x 200 ml
Ringer Solution	BT 5140.05	6 x 500 ml
RPMI 1640 + NaHCO ₃ + L-Glutamine + Phenol red	BT 5295.01	6 x 100 ml
RPMI 1640 + NaHCO ₃ + L-Glutamine + Phenol red	BT 5295.02	6 x 200 ml
RPMI 1640 + NaHCO ₃ + L-Glutamine + Phenol red	BT 5295.05	6 x 500 ml
Sabouraud Dextrose Broth	BT 5019.01	6 x 100 ml
Sabouraud Dextrose Broth	BT 5019.02	6 x 200 ml
Sabouraud Dextrose Broth	BT 5019.05	6 x 500 ml
Sabouraud Dextrose LAB-AGAR™	BT 5107.01	6 x 100 ml
Sabouraud Dextrose LAB-AGAR™	BT 5107.02	6 x 200 ml
Sabouraud Dextrose LAB-AGAR™	BT 5107.05	6 x 500 ml
Sabouraud Dextrose with Chloramphenicol LAB-AGAR™	BT 5108.01	6 x 100 ml
Sabouraud Dextrose with Chloramphenicol LAB-AGAR™	BT 5108.02	6 x 200 ml
Sabouraud Dextrose with Chloramphenicol LAB-AGAR™	BT 5108.05	6 x 500 ml
Sabouraud 2% Dextrose with Chloramphenicol LAB-AGAR™	BT 5109.01	6 x 100 ml
Sabouraud 2% Dextrose with Chloramphenicol LAB-AGAR™	BT 5109.02	6 x 200 ml
Sabouraud 2% Dextrose with Chloramphenicol LAB-AGAR™	BT 5109.05	6 x 500 ml
Sabouraud Dextrose with Chloramphenicol and Cycloheximide LAB-AGAR™	BT 5025.01	6 x 100 ml
Sabouraud Dextrose with Chloramphenicol and Cycloheximide LAB-AGAR™	BT 5025.02	6 x 200 ml
Sabouraud Dextrose with Chloramphenicol and Cycloheximide LAB-AGAR™	BT 5025.05	6 x 500 ml
Sabouraud Dextrose with Chloramphenicol and Gentamicin LAB-AGAR™	BT 5142.01	6 x 100 ml
Sabouraud Dextrose with Chloramphenicol and Gentamicin LAB-AGAR™	BT 5142.02	6 x 200 ml
Sabouraud Dextrose with Chloramphenicol and Gentamicin LAB-AGAR™	BT 5142.05	6 x 500 ml
Salicyne Broth	BT 5091.01	6 x 100 ml
Salicyne Broth	BT 5091.02	6 x 200 ml
Salicyne Broth	BT 5091.05	6 x 500 ml
Salmonella Shigella LAB-AGAR™	BT 5120.01	6 x 100 ml
Salmonella Shigella LAB-AGAR™	BT 5120.02	6 x 200 ml
Salmonella Shigella LAB-AGAR™	BT 5120.05	6 x 500 ml
Schaedler LAB-AGAR™ Base	BT 5013.01	6 x 100 ml
Schaedler LAB-AGAR™ Base	BT 5013.02	6 x 200 ml
Schaedler LAB-AGAR™ Base	BT 5013.05	6 x 500 ml
Schaedler Broth	BT 5014.01	6 x 100 ml
Schaedler Broth	BT 5014.02	6 x 200 ml
Schaedler Broth	BT 5014.05	6 x 500 ml
Schaedler Broth+ wit. K3	BT 5012.01	6 x 100 ml
Schaedler Broth+ wit. K3	BT 5012.02	6 x 200 ml
Schaedler Broth+ wit. K3	BT 5012.05	6 x 500 ml
Simmons Citrate LAB-AGAR™	BT 5101.01	6 x 100 ml
Simmons Citrate LAB-AGAR™	BT 5101.02	6 x 200 ml
Simmons Citrate LAB-AGAR™	BT 5101.05	6 x 500 ml
Sodium Selenite Broth	BT 5068.01	6 x 100 ml
Sodium Selenite Broth	BT 5068.02	6 x 200 ml
Sodium Selenite Broth	BT 5068.05	6 x 500 ml
Sodium Selenite Cystine Broth	BT 7013.01	6 x 100 ml
Sodium Selenite Cystine Broth	BT 7013.02	6 x 200 ml
Sodium Selenite Cystine Broth	BT 7013.05	6 x 500 ml
Sorbitol Broth	BT 5057.01	6 x 100 ml

Prepared culture media in bottles: 90, 100, 200, 225 and 500 ml

Description	Cat.No	Package
Sorbitol Broth	BT 5057.02	6 x 200 ml
Sorbitol Broth	BT 5057.05	6 x 500 ml
Sorbitol Mac Conkey LAB-AGAR™	BT 5021.01	6 x 100 ml
Sorbitol Mac Conkey LAB-AGAR™	BT 5021.02	6 x 200 ml
Sorbitol Mac Conkey LAB-AGAR™	BT 5021.05	6 x 500 ml
Thioglycollate Fluid Medium	BT 5128.01	6 x 100 ml
Thioglycollate Fluid Medium	BT 5128.02	6 x 200 ml
Thioglycollate Fluid Medium	BT 5128.05	6 x 500 ml
Todd Hewitt Broth	BT 5005.01	6 x 100 ml
Todd Hewitt Broth	BT 5005.02	6 x 200 ml
Todd Hewitt Broth	BT 5005.05	6 x 500 ml
Trehalose Broth	BT 5092.01	6 x 100 ml
Trehalose Broth	BT 5092.02	6 x 200 ml
Trehalose Broth	BT 5092.05	6 x 500 ml
Trehalose Broth (ISO 10273) (Yersinia)	BT 5241.01	6 x 100 ml
Trehalose Broth (ISO 10273) (Yersinia)	BT 5241.02	6 x 200 ml
Trehalose Broth (ISO 10273) (Yersinia)	BT 5241.05	6 x 500 ml
Trichomonas Broth	BT 5024.01	6 x 100 ml
Trichomonas Broth	BT 5024.02	6 x 200 ml
Trichomonas Broth	BT 5024.05	6 x 500 ml
Trypticasein Soy LAB-AGAR™	BT 5000.01	6 x 100 ml
Trypticasein Soy LAB-AGAR™	BT 5000.02	6 x 200 ml
Trypticasein Soy LAB-AGAR™	BT 5000.05	6 x 500 ml
Trypticasein Soy Broth	BT 5123.01	6 x 100 ml
Trypticasein Soy Broth	BT 5123.01	6 x 200 ml
Trypticasein Soy Broth	BT 5123.01	6 x 500 ml
Trypticasein Soy Broth	BT 5192	6 x 90 ml
Tryptophan Culture Broth	BT 5030.01	6 x 100 ml
Tryptophan Culture Broth	BT 5030.02	6 x 200 ml
Tryptophan Culture Broth	BT 5030.05	6 x 500 ml
TSYEA LAB – AGAR™ acc. to ISO 11290 (Listeria)	BT 5066.01	6 x 100 ml
TSYEA LAB – AGAR™ acc. to ISO 11290 (Listeria)	BT 5066.02	6 x 200 ml
TSYEA LAB – AGAR™ acc. to ISO 11290 (Listeria)	BT 5066.05	6 x 500 ml
TSYEB Broth	BT 5296.01	6 x 100 ml
TSYEB Broth	BT 5296.02	6 x 200 ml
TSYEB Broth	BT 5296.05	6 x 500 ml
Urea Broth Modified acc. to PN-A- 04023 (Hormeach and Munil)	BT 5040.01	6 x 100 ml
Urea Broth Modified acc. to PN-A- 04023 (Hormeach and Munil)	BT 5040.02	6 x 200 ml
Urea Broth Modified acc. to PN-A- 04023 (Hormeach and Munil)	BT 5040.05	6 x 500 ml
Urea Christensen Broth	BT 5028.01	6 x 100 ml
Urea Christensen Broth	BT 5028.02	6 x 200 ml
Urea Christensen Broth	BT 5028.05	6 x 500 ml
Urea Indol Broth	BT 5061.01	6 x 100 ml
Urea Indol Broth	BT 5061.02	6 x 200 ml
Urea Indol Broth	BT 5061.05	6 x 500 ml
VP Medium acc to. ISO 11290 (Listeria)	BT 5256.01	6 x 100 ml
VP Medium acc to. ISO 11290 (Listeria)	BT 5256.02	6 x 200 ml
VP Medium acc to. ISO 11290 (Listeria)	BT 5256.05	6 x 500 ml
XLD LAB-AGAR™	BT 5121.01	6 x 100 ml
XLD LAB-AGAR™	BT 5121.02	6 x 200 ml

Description	Cat.No	Package
XLD LAB-AGAR™	BT 5121.05	6 x 500 ml
Yersinia CIN LAB-AGAR™	BT 5035.01	6 x 100 ml
Yersinia CIN LAB-AGAR™	BT 5035.02	6 x 200 ml
Yersinia CIN LAB-AGAR™	BT 5035.05	6 x 500 ml

Prepared culture media in tubes**Milieux de culture prepares en tubes**

Description	Cat.No	Package
0,85 % Sodium Chloride	PW 3145	50 x 2 ml
0,85 % Sodium Chloride	PW 3053	50 x 5 ml
0,85 % Sodium Chloride	PW 3151	50 x 9 ml
0,85 % Sodium Chloride	PW 3232	50 x 10 ml
0,85 % Sodium Chloride	PW 3240	50 x 3 ml (Glass)
0,9 % Sodium Chloride	PW 3087	50 x 1 ml
0,9 % Sodium Chloride	PW 3150	50 x 2 ml
0,9 % Sodium Chloride	PW 3088	50 x 3 ml
0,9 % Sodium Chloride	PW 3054	50 x 5 ml
0,9 % Sodium Chloride	PW 3198	50 x 3 ml (Glass)
1% Glucose Broth	PW 3250	50 x 5 ml
5% Lactose Broth	PW 3132	50 x 5 ml
10% Lactose Broth	PW 3007	50 x 5 ml
10% Lactose Broth	PW 3099	50 x 3 ml
Acetamide Medium	PW 4214	50 x 5 ml
Adonitole Broth	PW 3223	50 x 5 ml
Alkaline Peptone Water	PW 3247	50 x 5 ml
Alkaline Peptone Water(10x)	PW 3248	50 x 5 ml
Arabinose Broth	PW 3046	50 x 5 ml
Arabinose Broth	PW 3169	50 x 7 ml
Arginine acc. to Falkow Broth	PW 3048	50 x 5 ml
Arginine acc. to Falkow Broth	PW 3184	50 x 3 ml
Azide Dextrose Broth	PW 3162	50 x 5 ml
Bile Esculin LAB-AGAR™	PW 3104	50 x 3 ml (Slant)
Bile salts	PW 3249	50 x 5 ml
Brain Heart Infusion Broth	PW 3015	50 x 10 ml
Brain Heart Infusion Broth	PW 3094	50 x 9 ml
Brain Heart Infusion Broth	PW 3031	50 x 5 ml
Brain Heart Infusion Broth	PW 4202	50 x 3 ml
Brain Heart Infusion Broth	PW 3154	50 x 4 ml
Brain Heart Infusion Broth	PW 3155	50 x 9 ml (Glass)
Brain Heart Infusion Broth	PW 3174	50 x 10 ml (Glass)
Brain Heart Infusion Broth	PW 3196	50 x 5 ml (Glass)
Brucella Broth	PW 3097	50 x 5 ml
Brucella Broth	PW 3172	50 x 10 ml (Glass)
Brucella Broth	PW 3200	50 x 3 ml (Glass)
Brucella Broth + wit. K3	PW 3148	50 x 5 ml
Bulion (Broth)+ 6,5% NaCl	PW 3039	50 x 5 ml
Bulion (Broth) pH 9,6	PW 3076	50 x 5 ml
Bulion zwykły (Nutritious Broth)	PW 3037	50 x 5 ml
Bulion zwykły (Nutritious Broth)	PW 3038	50 x 10 ml
Celobiose Broth	PW 3080	50 x 5 ml

Prepared culture media in tubes
Milieux de culture prepares en tubes

Description	Cat.No	Package
Chocolate LAB-AGAR™ + PV	PW 3251	50 x 5 ml (Butt)
Chromogenic Malassezia LAB-AGAR™	PW 3252	50 x 5 ml (Butt)
Chromogenic Strepto B Lab-AGAR™	PW 3239	50 x 5 ml (Butt)
Citrate Christensen LAB-AGAR™	PW 3056	50 x 5 ml (Slant)
Dermatophytes LAB-AGAR™	PW 3059	50 x 5 ml (Slant)
Distilled Water	PW 3089	50 x 1 ml
Distilled Water	PW 3090	50 x 3 ml
Distilled Water	PW 3091	50 x 5 ml
Distilled Water	PW 3160	50 x 10 ml
Dulcitol Broth	PW 3049	50 x 5 ml
Dulcitol Broth	PW 3165	50 x 7 ml
EMB LAB-AGAR™	PW 3253	50 x 5 ml (Butt)
Enrichment Broth	PW 3254	50 x 5 ml
Enrichment LAB-AGAR™	PW 3255	50 x 5 ml (Slant)
Enterococcosel LAB-AGAR™	PW 3256	50 x 5 ml (Butt)
Enterococcus Arabinose Broth	PW 3063	50 x 5 ml
Ewing Malonate Modified Broth	PW 3020	50 x 5 ml
Ewing Malonate Modified Broth	PW 3167	50 x 3 ml
Ewing Malonate Modified Broth	PW 3175	50 x 2 ml
Eugon Broth	PW 3257	50 x 5 ml
Falkow Medium Base	PW 3010	50 x 5 ml
Falkow Medium Base	PW 3103	50 x 3 ml
Fungisel LAB-AGAR™	PW 3258	50 x 5 ml(Slant)
Glicerol Broth	PW 3081	50 x 5 ml
Glucose Enrichment Broth	PW 3025	50 x 5 ml
Glucose Enrichment Broth	PW 3092	50 x 9 ml
Hugh Leifson Medium	PW 3105	50 x 5 ml
Inositol Broth	PW 3082	50 x 5 ml
KCN Broth Base	PW 3078	50 x 5 ml
King B LAB-AGAR™	PW 3109	50 x 7 ml (Slant)
Kligler Iron LAB-AGAR™	PW 3016	50 x 5 ml (Butt)
Kligler Iron LAB-AGAR™	PW 3057	50 x 5 ml(Slant)
Kligler Iron LAB-AGAR™	PW 3116	50 x 3 ml (Slant)
Ksylose Broth	PW 3021	50 x 5 ml
Ksylose Broth acc to ISO 11290 (Listeria)	PW 3000	50 x 5 ml
Liver Broth	PW 3122	50x10 ml
Lysine acc.to Falkow Broth	PW 3026	50 x 5 ml
Lysine acc.to Falkow Broth	PW 3118	50 x 9 ml
Lysine acc.to Falkow Broth	PW 3102	50 x 3 ml
Lysine Iron LAB-AGAR™	PW 3195	50 x 7 ml(Slant)
Maltose Broth	PW 3075	50 x 5 ml
Mannitol Broth	PW 3176	50 x 2 ml
Mannitol Broth	PW 3044	50 x 5 ml
Mannitol Broth	PW 3163	50 x 7 ml
Mannitol Salt LAB-AGAR™	PW 3259	50 x 5 ml (Butt)
Melibiose Broth (ISO 10273) (Yersinia)	PW 3231	50 x 5 ml
Melibiose Broth	PW 3086	50 x 5 ml
MIU Medium	PW 3191	50 x 3 ml
Modified Motility LAB-AGAR™	PW 3096	50 x 5 ml
Modified Motility LAB-AGAR™	PW 3106	50 x 3 ml

Prepared culture media in tubes

Milieux de culture prepares en tubes

Description	Cat.No	Package
Motility LAB-AGAR™	PW 3023	50 x 5 ml
Motility LAB-AGAR™ acc. to ISO 11290 (Listeria)	PW 3221	50 x 5 ml
MRVP Medium	PW 3036	50 x 5 ml
Mueller Hinton Broth	PW 3028	50 x 5 ml
Mueller Hinton Broth	PW 3142	50 x 6 ml
Mueller Hinton II Broth	PW 3041	50 x 11 ml
Mueller Hinton Broth with Laked Horse Blood	PW 3051	50 x 11 ml
Nitrate Broth	PW 3047	50 x 5 ml
Nutrient Broth	PW 3194	50 x 9 ml
Nutrient Broth	PW8001	50 x 5 ml
Nutrient No 2 Broth	PW 3003	50 x 3 ml
Nutrient No 2 Broth	PW 3004	50 x 5 ml
Nutrient No 2 Broth	PW 3079	50 x 7 ml
Nutrient No 2 Broth	PW 3245	50 x 9 ml
Nutrient No 2 Broth	PW 3246	50x10 ml
Nutrient LAB-AGAR™	PW 3065	50 x 5 ml (Slant)
Nutrient LAB-AGAR™ at pH 7,0	PW 3260	50 x 5 ml (Slant)
Nutrient No 1 LAB-AGAR™	PW 3179	50 x 7 ml (Slant)
Nutrient No 2 LAB-AGAR™	PW 3111	50 x 5 ml (Slant)
Nutrient No 2 LAB-AGAR™	PW 3218	50 x 9 ml (Slant)
Nutrient Pneumococcus Broth	PW 3029	50 x 1 ml
Nutrient S LAB-AGAR™ (for bacterial strains vitality preservation)	PW 3134	50 x 3 ml (Butt)
Nutrient S LAB-AGAR™ (for bacterial strains vitality preservation)	PW 3135	50 x 1 ml (Butt)
Nutrient S LAB-AGAR™ (for bacterial strains vitality preservation)	PW 3220	50 x 5 ml (butt)
Ornithine acc. to Falkow Broth	PW 3009	50 x 5 ml
Ornithine acc. to Falkow Broth	PW 3168	50 x 7 ml
Ornithine acc. to Falkow Broth	PW 3177	50 x 2 ml
Ornithine acc. to Falkow Broth	PW 3185	50 x 3 ml
Peptone Water	PW 3074	50 x 5 ml
Phenylalanine LAB-AGAR™	PW 3024	50 x 5 ml (Butt)
Phenylalanine LAB-AGAR™	PW 3061	50 x 5 ml (Slant)
Phenylalanine LAB-AGAR™	PW 3119	50 x 3 ml (Slant)
Potato Dextrose LAB-AGAR™	PW 3261	50 x 5 ml (Butt)
Rafinose Broth	PW 3083	50 x 5 ml
RambaQUICK LIM STREP B Broth	PW 3262	50 x 5 ml
RambaQUICK Salmonella Broth	PW 3263	50 x 5 ml
Ramnose Broth acc. to ISO 11290 (Listeria)	PW 3001	50 x 5 ml
Ramnose w. Andrade Indicator Broth	PW 3224	50 x 5 ml
Ringer Solution	PW 3107	50 x 9 ml
Ringer Solution	PW 3107S	50 x 9 ml (Glass)
Ringer Solution	PW 3164	50 x 1,5 ml
Sabouraud Dextrose Broth	PW 3022	50 x 10 ml
Sabouraud Dextrose Broth	PW 3034	50 x 5 ml
Sabouraud Dextrose Broth	PW 3207	50 x 1 ml
Sabouraud Dextrose Broth	PW 3181	50 x 10 ml (Glass)
Sabouraud with Chloramphenicol Broth	PW 3241	50 x 9 ml
Sabouraud Dextrose Broth + CH + G	PW 3264	50 x 5 ml
Sabouraud Dextrose with Chloramphenicol LAB-AGAR™	PW 3112	50 x 10ml (Slant)
Sabouraud Dextrose with Chloramphenicol and Cycloheximide LAB-AGAR™	PW 3064	50 x 5ml (Slant)
Sabouraud Dextrose with Chloramphenicol and Gentamicin LAB-AGAR™	PW 3206	50 x 5ml (Slant)

Prepared culture media in tubes
Milieux de culture prepares en tubes

Description	Cat.No	Package
Sacharose Broth	PW 3225	50 x 5 ml
Salicyne Broth	PW 3084	50 x 5 ml
Schaedler Broth	PW 3017	50 x 10 ml
Schaedler Broth	PW 3069	50 x 5 ml
Schaedler Broth	PW 3173	50 x 10 ml (Glass)
Schaedler Broth+ wit K3	PW 3018	50 x 10 ml
Schaedler Broth +wit. K3	PW 3014	50 x 5 ml
Schaedler Broth +wit. K3	PW 3035	50 x 8 ml
Schaedler Broth +wit. K3	PW 3201	50 x 13 ml
Schaedler Modified Broth	PW 3090	50 x 7 ml
Schaedler Modified Broth	PW 3204	50 x 10 ml
Simmons Citrate LAB-AGAR™	PW 3019	50 x 5 ml (Butt)
Simmons Citrate LAB-AGAR™	PW 3058	50 x 5ml (Slant)
Simmons Citrate LAB-AGAR™	PW 3117	50 x 3ml (Slant)
Slanetz and Bartley LAB-AGAR™	PW 3265	50 x 5 ml (Butt)
Sodium Cystine Broth	PW 3192	50 x 9 ml
Sodium Selenite Broth	PW 3189	50 x 3 ml
Sodium Selenite Broth	PW 4209	50 x 5 ml
Sodium Selenite Broth	PW 3040R	50 x 6 ml
Sodium Selenite Broth	PW 3040	50 x 6 ml
Sodium Selenite Broth (transport tube)	PW 3040.T	50 x 10 ml
Sodium Selenite Broth	PW 3030	50 x 7 ml
Sodium Selenite Broth	PW 3095	50 x 9 ml
Sodium Selenite Broth	PW 3008	50 x 10 ml
Sorbitol Broth	PW 3045	50 x 5 ml
Sorbitol Broth	PW 3166	50 x 7 ml
Thioglycollate Fluid Medium	PW 3033	50 x 5 ml
Thioglycollate Fluid Medium	PW 4004	50 x 10 ml
Thioglycollate Fluid Medium	PW 3153	50 x 9 ml
Thioglycollate Fluid Medium	PW 3157	50 x 9 ml (Glass)
Thioglycollate Fluid Medium	PW 3187	50 x 10 ml (Glass)
Todd Hewitt Broth	PW 3027	50 x 5 ml
Todd Hewitt Broth	PW 3077	50 x 10 ml
Todd Hewitt Broth	PW 3188	50 x 3 ml
Todd Hewitt Broth with Nalidixic Acid and Colistin	PW 3203	50 x 5 ml
Todd Hewitt Broth with Nalidixic Acid and Gentamicin	PW 3050	50 x 5 ml
Todd Hewitt Broth with Nalidixic Acid and Gentamicin	PW 3171	50 x 2 ml
Trehalose Broth	PW 3085	50 x 5 ml
Trehalose Broth (ISO 10273) (Yersinia)	PW 3230	50 x 5 ml
Trichomonas Broth	PW 3011	50 x 3 ml
Trichomonas Broth	PW 3180	50 x 2,5 ml
Triple Sugar Iron LAB-AGAR™	PW 4204	50 x 7ml (Butt)
Triple Sugar Iron LAB-AGAR™	PW 3060	50 x 7ml (Slant)
Triple Sugar Iron LAB-AGAR™	PW 3193	50 x 9ml (Slant)
Trypticasein Soy Broth	PW 3101	50 x 3 ml
Trypticasein Soy Broth	PW 4019	50 x 10 ml
Trypticasein Soy Broth	PW 4000	50 x 5 ml
Trypticasein Soy Broth	PW 3152	50 x 9 ml
Trypticasein Soy Broth	PW 3156	50 x 9 ml (Glass)
Trypticasein Soy Broth	PW 3197	50 x 5 ml (Glass)

Prepared culture media in tubes

Milieux de culture préparés en tubes

Description	Cat.No	Package
Tryptophan Culture Broth	PW 3178	50 x 2 ml
Tryptophan Culture Broth	PW 3098	50 x 3 ml
Tryptophan Culture Broth	PW 3012	50 x 5 ml
TSA LAB-AGAR™	PW 3062	50 x 5 ml (Slant)
TSA LAB-AGAR™	PW 3146	50 x 10 ml (Slant)
TSA LAB-AGAR™	PW 3205	50 x 7 ml (Slant)
TSYEB Broth (Listeria)	PW 3147	50 x 5 ml
Urea Broth Modified acc. to PN-A- 04023 (Hormeach and Munil)	PW 3100	50 x 3 ml
Urea Broth Modified acc. to PN-A- 04023 (Hormeach and Munil)	PW 3002	50 x 5 ml
Urea Chritensen Broth	PW 3006	50 x 5 ml
Urea Chritensen Broth	PW 3227	50 x 3 ml
Urea Indol Broth	PW 4212	50 x 5 ml
Urea Indol Broth	PW 3149	50 x 3 ml
Urea Indol Broth	PW 3242	50 x 9 ml
Urea LAB-AGAR™	PW 4233	50 x 5 ml (Slant)
Urea LAB-AGAR™	PW 3120	50 x 7 ml (Slant)
Urea LAB-AGAR™	PW 3183	50 x 10 ml (Slant)
VP Medium acc to. ISO 11290 (Listeria)	PW 3236	50 x 3 ml

Media in tubes for diagnosis of mycobacteria

Milieux de culture en tubes pour diagnostic des mycobactéries

Description	Cat.No	Package
Lowenstein - Jensen Medium w/o MG	153/T	20 tubes*
Lowenstein - Jensen Medium + G + PACT	152/T	20 tubes*
Lowenstein - Jensen Medium + Ofloxacin 5mg/l	151/T	20 tubes*
Lowenstein - Jensen Medium + Ofloxacin 2,5mg/l	150/T	20 tubes*
Lowenstein - Jensen Medium + Rifampicin 80mg/l	149/T	20 tubes*
Lowenstein - Jensen Medium + Rifampicin 40mg/l	148/T	20 tubes*
Lowenstein - Jensen Medium + Ethambutol 4mg/l	147/T	20 tubes*
Lowenstein - Jensen Medium + Ethambutol 2mg/l	146/T	20 tubes*
Lowenstein - Jensen Medium + Isoniazid 0,4mg/l	145/T	20 tubes*
Lowenstein - Jensen Medium + Isoniazid 0,2mg/l	144/T	20 tubes*
Lowenstein - Jensen Medium + Streptomycin 8mg/l	143/T	20 tubes*
Lowenstein - Jensen Medium + Streptomycin 4mg/l	142/T	20 tubes*
Lowenstein - Jensen Medium + MG without Glycerol	141/T	20 tubes*
Lowenstein - Jensen Medium + MG + Glycerol	140/T	20 tubes*
Lowenstein - Jensen Medium + MG + Glycerol + P	139/T	20 tubes*
Lowenstein - Jensen Medium + MG + Glycerol + TCH	138/T	20 tubes*
Stonebrink Medium + MG + P	137/T	20 tubes*
Stonebrink Medium + MG + PACT	135/T	20 tubes*

*available in glass tubes: 5 ml, 7 ml, 9 ml, 10 ml

ABV- Amphotericin+Bacitracin+Vancomycin

HB- blood horse

SB - blood sheep

C- Chloramphenicol

CH- Chloramphenicol

CN - Cetrimide+Nalidixid acid

CIN- Cefsulodin+Irgasan+Novobiocin

CV- crystal violet

G - Gentamycin

Ge - Gentamycin

MG-Malachite Green

NEO- Neomycin P- Pyruvate

PACT- Polymyxin+Amphotericin+Carbenicillin+Trimethoprim

PV- PolyVitex

TCH- Thiocarbohydrazide

V- Vancomycin

VA- Vancomycin

VCAT- Vancomycin+Colistin+Amphotericin+Trimethoprim

VCN - Vancomycin+Colistin+Nystatin

VCNT - Vancomycin+Colistin+Nystatin+Trimethoprim

Dehydrated culture media

Millieux de culture deshydrates

Description	Cat.No	Package
10 % Peptone Water	PS 524	500 g
10% Lactose Broth	PS 503	500 g
Acetamide Broth - Twin Pack	PS 601E	500 g
Acetamide Medium (bulion z acetamidem)	PS 214	500 g
Acetamide Twin Pack LAB-AGAR™	PS 600E	500 g
Acetate Differential LAB-AGAR	PS 263	500 g
Actinomyces Isolation LAB-AGAR™	PS 585	500 g
Aeromonas (Ryan) Medium Base	PS 602E	500 g
Aeromonas LAB AGAR™	PS 508	500 g
Alicyclobacillus Medium	PS 603E	500 g
Alkaline Peptone Water	PS 604E	500 g
Alkaline Peptone Water	PS 151	500 g
Alkaline Peptone Water (ASPW)	PS 282	500 g
Amies Transport LAB-AGAR™	PS 46	500 g
Amies w. Charcoal Transport Medium	PS 606E	500 g
Anaerobic (Brewer) LAB-AGAR™	PS 607E	500 g
Andrade Lactose Medium	PS 608E	500 g
Antibiotic Assay Medium No 1	PS 875E	500 g
Antibiotic Assay Medium No 10	PS 879E	500 g
Antibiotic Assay Medium No 11	PS 877E	500 g
Antibiotic Assay Medium No 13	PS 874E	500 g
Antibiotic Assay Medium No 19	PS 878E	500 g
Antibiotic Assay Medium No 2	PS 872E	500 g
Antibiotic Assay Medium No 3	PS 262	500 g
Antibiotic Assay Medium No 3	PS 881E	500 g
Antibiotic Assay Medium No 32	PS 873E	500 g
Antibiotic Assay Medium No 37	PS 884E	500 g
Antibiotic Assay Medium No 41	PS 880E	500 g
Antibiotic Assay Medium No 5	PS 876E	500 g
Antibiotic Medium No 1	PS 510	500 g
Antibiotic Medium No 2	PS 559	500 g
AOAC Lethen Broth	PS 217	500 g
APT Broth	PS 610E	500 g
APT LAB-AGAR™	PS 609E	500 g
Arginine Dihydrolase	PS 1038E	500 g
Asparagine Proline Broth	PS 611E	500 g
Aspergillus Differentiation Medium Base	PS 612E	500 g
Assay at pH 6,0 LAB-AGAR™	PS 206	500 g
Assay at pH 7,2 LAB-AGAR™	PS 207	500 g
Assay at pH 8,0 LAB-AGAR™	PS 208	500 g
Azide Blood LAB-AGAR™ Base	PS 613E	500 g
Azide Blood LAB-AGAR™ Base	PS 122	500 g
Azide Dextrose Rothe Broth	PS 614E	500 g
Azide Dextrose w.BCP Rothe Broth	PS 615E	500 g
B.C.P. Glucose LAB-AGAR™	PS 154	500 g
B.G.A. LAB-AGAR™ acc. to ISO 6579	PS 98	500 g
Bacillus cereus acc. to Mossel LAB-AGAR™ Base	PS 616E	500 g
Bacillus cereus Selective LAB-AGAR™ Base	PS 47	500 g

Dehydrated culture media

Millieux de culture deshydrates

Description	Cat.No	Package
Baird Parker LAB-AGAR™ Base	PS 617E	500 g
Baird Parker LAB-AGAR™ Base acc.to ISO 6888-1	PS 33	500 g
Baird Parker LAB-AGAR™ Base acc.to ISO 6888-2	PS 164	500 g
BAT Broth	PS 298	500 g
BAT LAB-AGAR™	PS 289	500 g
BIGGY LAB-AGAR™	PS 591	500 g
BIGGY LAB-AGAR™	PS 618E	500 g
Bile Esculin Azide LAB-AGAR™	PS 620E	500 g
Bile Esculin Azide LAB-AGAR™ acc.to ISO 7899-2	PS 88	500 g
Bile Esculin LAB-AGAR™ acc. to ISO 10773	PS 619E	500 g
Bile Esculin LAB-AGAR™	PS 155	500 g
Bismuth Sulfite LAB-AGAR™	PS 01	500 g
Bismuth Sulphite LAB-AGAR™ acc. to USP	PS 621E	500 g
Bismuth Sulphite LAB-AGAR™ Modified	PS 622E	500 g
Blickfeldta LAB-Agar Base	PS 348	500 g
Blood Free Campylobacter LAB-AGAR™ Base	PS 625E	500 g
Blood Free Campylobacter CCDA LAB-AGAR™ Base	PS 159	500 g
Blood LAB-AGAR™ Base	PS 02	500 g
Blood No 2 LAB-AGAR™ Base	PS 624E	500 g
Blood No 2 LAB-AGAR™ Base	PS 205	500 g
Bolton Broth Base	PS 626E	500 g
Bolton Broth Base	PS 179	500 g
Bordet Gengou LAB-AGAR™ Base	PS 627E	500 g
BPL LAB-AGAR™	PS 628E	500 g
Brain Heart Infusion 1 % LAB-AGAR™	PS 630E	500 g
Brain Heart Infusion Broth	PS 04	500 g
Brain Heart Infusion Broth	PS 631E	500 g
Brain Heart Infusion LAB-AGAR™	PS 629E	500 g
Brain Heart Infusion LAB-AGAR™	PS 03	500 g
Brilliant Green acc.to Kauffman LAB-AGAR™	PS 48	500 g
Brilliant Green Bile Broth 2%	PS 27	500 g
Brilliant Green Bile Broth 2%	PS 635E	500 g
Brilliant Green Bile Broth 2% acc. to ISO 4831	PS 636E	500 g
Brilliant Green LAB-AGAR™ acc. to ISO 6340	PS 634E	500 g
Brilliant Green Lab-AGAR™ Modified	PS 632E	500 g
Brilliant Green Sulpha LAB-AGAR™ Agar	PS 1036E	500 g
BRILLIANT GREEN TETRATHIONATE BILE BROTH (European Pharm.)	PS 199	500 g
Brilliant Green w. Phosphates LAB-AGAR™ Base	PS 633E	500 g
BRILLIANT GREEN XYLOSE LAB-AGAR™	PS 237	500 g
Brilliant Green, Phenol Red, Lactose Monohydrate, Sucrose LAB-AGAR™	PS 1037E	500 g
Bromocresol Purple Azide Broth	PS 637E	500 g
Bromocresol Purple Broth Base	PS 212	500 g
Bromocresol Purple Broth Base (ISO 11290)	PS 537	500 g
Bromocresol Purple w. Lactose Broth	PS 638E	500 g
Brucella Broth	PS 137	500 g
Brucella Broth Base	PS 640E	500 g
Brucella Lab-AGAR™ Base	PS 639E	500 g
Brucella LAB-AGAR™	PS 05	500 g
Bryant-Burkey Broth Base	PS 114	500 g
Buffered Charcoal Yeast Extract LAB-AGAR™ Base	PS 641E	500 g

Dehydrated culture media
Millieux de culture deshydrates

Description	Cat.No	Package
Buffered Peptone Water	PS 52	500 g
Buffered Peptone Water acc. to Eur. Pharm.	PS 181	500 g
Buffered Saline Peptone Water	PS 149	500 g
Buffered Yeast Extract LAB-AGAR™	PS 644E	500 g
Burkholderia cepacia LAB-AGAR™ Base	PS 645E	500 g
Burkholderia Cepacia LAB-AGAR™ Base	PS 254	500 g
Cameronia LAB-AGAR™	PS 246	500 g
Campylobacter acc. to Preston LAB-AGAR™ Base	PS 160	500 g
Campylobacter Enrichment Broth Base	PS 649E	500 g
Campylobacter LAB-AGAR™ Base	PS 648E	500 g
Carbohydrate Fermentation Broth Base	PS 650E	500 g
Cary-Blair Medium	PS 249	500 g
Cary-Blair Medium Base	PS 651E	500 g
Casman Broth Base	PS 653E	500 g
Casman LAB-AGAR™	PS 652E	500 g
Cetrimide Broth Base	PS 655E	500 g
Cetrimide LAB-AGAR™ Base	PS 654E	500 g
Cetrimide LAB-AGAR™	PS 49	500 g
Charcoal LAB-AGAR™ Base	PS 174	500 g
Charcoal LAB-AGAR™ Base	PS 656E	500 g
Charcoal with Niacin LAB-AGAR™	PS 657E	500 g
China Blue Lactose LAB-AGAR™	PS 658E	500 g
Chloramphenicol LAB-AGAR™	PS 106	500 g
Chloramphenicol Yeast Glucose LAB-AGAR™	PS 659E	500 g
CHO Medium Base	PS 660E	500 g
Cholera Medium Base	PS 661E	500 g
Christensen Citrate LAB-AGAR™	PS 662E	500 g
Chromogenic aureus LAB-AGAR™ Base	PS 663E	500 g
Chromogenic Bacillus LAB-AGAR™ Base	PS 664E	500 g
Chromogenic Candida LAB-AGAR™	PS 665E	500 g
Chromogenic Candida LAB-AGAR™	PS 590	500 g
Chromogenic Coliform (CCA) LAB-AGAR™ acc. to ISO 9308-1	PS 666E	500 g
Chromogenic Coliform LAB-AGAR acc. to ISO 9308-1	PS 267	500 g
Chromogenic Coliform w. SLS LAB-AGAR™	PS 667E	500 g
Chromogenic Cronobacter Isolation (CCI) LAB-AGAR™	PS 284	500 g
Chromogenic E. coli (TBX) LAB-AGAR™	PS 668E	500 g
Chromogenic E. coli LAB-AGAR™	PS 669E	500 g
Chromogenic E.coli O157 LAB-AGAR™	PS 589	500 g
Chromogenic ECC LAB-AGAR™	PS 670E	500 g
Chromogenic Enterobacter sakazakki LAB-AGAR™	PS 671E	500 g
Chromogenic Enterococcus faecium LAB-AGAR™ Base	PS 673E	500 g
Chromogenic Enterococcus LAB-AGAR™	PS 672E	500 g
Chromogenic Klebsiella LAB-AGAR™ Base	PS 674E	500 g
Chromogenic Listeria acc.to Ottaviani and Agostii LAB-AGAR™ Base	PS 165	500 g
Chromogenic Listeria LAB-AGAR™ Base	PS 605E	500 g
Chromogenic Listeria LAB-AGAR™ Base Modified	PS 675E	500 g
Chromogenic Pseudomonas LAB-AGAR™	PS 224	500 g
Chromogenic Salmonella LAB-AGAR™	PS 598	500 g
Chromogenic Salmonella LAB-AGAR™	PS 676E	500 g
Chromogenic UTI LAB-AGAR™	PS 677E	500 g

Dehydrated culture media

Milieux de culture deshydrates

Description	Cat.No	Package
Chromogenic Vibrio LAB-AGAR™	PS 274	500 g
Citrate acc.to Christensen LAB AGAR™	PS 530	500 g
Clark Medium	PS 253	500 g
Clausen Medium	PS 678E	500 g
CLED LAB-AGAR™	PS 34	500 g
CLED LAB-AGAR™	PS 646E	500 g
CLED LAB-AGAR™ w. Andrade's Indicator	PS 647E	500 g
CLED with Andrade Indicator LAB-AGAR™	PS 113	500 g
Clostridium Broth	PS 1039E	500 g
Clostridium Difficile LAB-AGAR™ Base	PS 679E	500 g
Clostridium difficile LAB-AGAR™ Base	PS 176	500 g
Columbia CNA LAB-AGAR™	PS 190	500 g
Columbia LAB-AGAR™ Base	PS 06	500 g
Cooked Meat Medium	PS 682E	500 g
Corn Meal LAB-AGAR™	PS 683E	500 g
Cronobacter Selective Broth (CSB) Base	PS 276	500 g
Cystine Heart LAB-AGAR™ Base	PS 684E	500 g
Czapek Dox Broth	PS 686E	500 g
Czapek Dox LAB-AGAR™	PS 685E	500 g
Czapek Dox LAB-AGAR™	PS 28	500 g
Czapek-Dox Modified Broth	PS 280	500 g
D.C.L.S. LAB-AGAR™	PS 687E	500 g
D.T.M. LAB-AGAR™ Base	PS 688E	500 g
D/E Neutralising Broth (acc.to USP)	PS 187	500 g
Decarboxylase Broth Base Moeller	PS 689E	500 g
Decarboxylase Test Medium Base Falkow	PS 690E	500 g
Dermatophytes LAB AGAR™	PS 586	500 g
Desoxycholate Citrate LAB-AGAR™	PS 693E	500 g
Desoxycholate Citrate LAB-AGAR™ Modified	PS 694E	500 g
Desoxycholate Citrate LAB-AGAR™	PS 39	500 g
Desoxycholate LAB-AGAR™	PS 691E	500 g
Desoxycholate Lactose LAB-AGAR™	PS 695E	500 g
Desoxycholate Lactose LAB-AGAR™	PS 40	500 g
Desoxycholate Leifson LAB-AGAR™ Modified	PS 692E	500 g
DEV Lactose Broth	PS 696E	500 g
Dextrose Broth	PS 697E	500 g
Dextrose Peptone Broth	PS 699E	500 g
Dextrose Peptone LAB-AGAR™	PS 698E	500 g
Dextrose Proteose Peptone LAB-AGAR™ Base	PS 700E	500 g
Dextrose Tryptone Broth	PS 702E	500 g
Dextrose Tryptone LAB-AGAR™	PS 701E	500 g
Diagnostic Sensitivity Test LAB-AGAR™	PS 703E	500 g
Diagnostic Sensitivity Test LAB-AGAR™	PS 230	500 g
Dichloran Glycerol DG 18 LAB-AGAR™	PS 168	500 g
Dichloran Glycerol Medium Base	PS 704E	500 g
Dichloran Rose Bengal Chloramphenicol (DRBC) LAB-AGAR™	PS 706E	500 g
Dichloran w. Rose Bengal Medium Base	PS 705E	500 g
Differential Clostridial Broth Modified	PS 707E	500 g
Differential Reinforced Clostridial Broth Base	PS 708E	500 g
DNase Test LAB-AGAR™ Base	PS 709E	500 g

Dehydrated culture media
Millieux de culture deshydrates

Description	Cat.No	Package
DNAse Test LAB-AGAR™	PS 50	500 g
DRBC LAB-AGAR™	PS 136	500 g
DRCM Broth (Differential Clostridial Medium)	PS 220	500 g
Drigalski Li™us Lactose LAB-AGAR™	PS 710E	500 g
E. coli Chromogenic Medium	PS 121	500 g
E.S.T.Y Medium	PS 80	500 g
EC Broth	PS 711E	500 g
EC Medium	PS 51	500 g
Edwards LAB-AGAR™	PS 552	500 g
Edward's Medium Modified	PS 712E	500 g
EE Broth	PS 171	500 g
EE Broth Mossel	PS 713E	500 g
Eijkman Lactose Broth	PS 714E	500 g
Eijkman Lactose Medium	PS 76	500 g
EMB LAB-AGAR™	PS 715E	500 g
EMB Levin LAB-AGAR™	PS 716E	500 g
Endo LAB-AGAR™	PS 717E	500 g
Endo LAB-AGAR™ Base	PS 718E	500 g
Endo LAB-AGAR™	PS 29	500 g
Endo LES LAB-AGAR™ Base	PS 89	500 g
Endo LES with Fuchsine LAB-AGAR™	PS 195	500 g
Enrichment Broth	PS 110	500 g
Enrichment LAB-AGAR™	PS 109	500 g
Enterococcus Confirmatory LAB-AGAR™	PS 95	500 g
Eosin Methylene Blue LAB-AGAR™	PS 41	500 g
ESIA LAB-AGAR™ (Enterobacter sakazakii)	PS 173	500 g
Ethyl Violet Azide EVA Broth	PS 719E	500 g
EUGON LT 100 Broth Base	PS 293	500 g
Eugon LAB-AGAR™	PS 42	500 g
Eugonic Broth	PS 721E	500 g
Eugonic LAB-AGAR™	PS 720E	500 g
Eugonic LT 100 Broth Base	PS 722E	500 g
EVA Broth	PS 83	500 g
Ewing Malonate Modified Broth	PS 53	500 g
Fecal Coliforms Broth Base	PS 93	500 g
Fecal Coliforms LAB-AGAR™ Base	PS 87	500 g
Folic Acid Assay Medium	PS 723E	500 g
Fraser Broth Base	PS 101	500 g
Fraser Broth Base	PS 725E	500 g
Fraser Listeria Enrichment Broth Base	PS 99	500 g
Fraser Secondary Enrichment Broth Base	PS 726E	500 g
G.C. LAB-AGAR™ Base	PS 727E	500 g
G.C. LAB-AGAR™ Base	PS 139	500 g
Gassner Lactose LAB-AGAR™	PS 728E	500 g
GBS Medium Base	PS 729E	500 g
Gelatin LAB-AGAR™ Agar	PS 1034E	500 g
Gelatin Lactose Medium	PS 112	500 g
Giolitti-Cantoni Broth Base with Polysorbate	PS 731E	500 g
Giolitti-Cantoni with Polysorbate Broth	PS 148	500 g
Giolitti-Cantoni without Polysorbate Broth	PS 54	500 g

Dehydrated culture media

Millieux de culture deshydrates

Description	Cat.No	Package
Glucose Broth	PS 733E	500 g
Glucose Chloramphenicol Broth	PS 117	500 g
Glucose LAB-AGAR™	PS 732E	500 g
Glucose OF Medium	PS 268	500 g
GN Hajna Broth	PS 734E	500 g
GSP LAB-AGAR™	PS 223	500 g
Haemophilus Test LAB-AGAR™ Base	PS 735E	500 g
Haemophilus Test LAB-AGAR™ Base	PS 182	500 g
Half Fraser Broth Base	PS 100	500 g
Hartley's Digest Broth	PS 736E	500 g
HC LAB-AGAR™ BASE	PS 216	500 g
Heart Infusion LAB-AGAR™	PS 737E	500 g
Hektoen Enteric LAB-AGAR™	PS 738E	500 g
Hektoen LAB-AGAR™	PS 07	500 g
High Salt Nutrient LAB-AGAR™	PS 739E	500 g
High Salt Peptone Yeast Extract LAB-AGAR™	PS 740E	500 g
Hoyle's Medium Base	PS 743E	500 g
Hugh Leifson w. Glucose Medium	PS 744E	500 g
Indol Nitrate Medium	PS 55	500 g
Indol Nitrate Medium	PS 745E	500 g
Inositol Brilliant Green Bile LAB-AGAR™	PS 746E	500 g
Iron Sulphite LAB-AGAR™	PS 747E	500 g
Islams GBS Medium Base	PS 748E	500 g
ISO-SENSITEST LAB-AGAR™	PS 239	500 g
ITC Broth Base	PS 749E	500 g
ITC Broth Base	PS 169	500 g
IUT Mycobacteria Medium Base	PS 750E	500 g
K.R.A.N.E.P. LAB-AGAR™ Base	PS 751E	500 g
KAA Confirmatory LAB-AGAR™	PS 123	500 g
Kanamycin Esculin Azide Broth	PS 753E	500 g
Kanamycin Esculin Azide LAB-AGAR™	PS 752E	500 g
Karmali Campylobacter LAB-AGAR™ Base	PS 754E	500 g
Karmali LAB-AGAR™ Base	PS 178	500 g
KCN Broth Base	PS 545	500 g
KF Streptococcal acc. to Kranter Broth Base	PS 756E	500 g
KF Streptococcal acc. to Kranter LAB-AGAR™ Base	PS 755E	500 g
KF Streptococcal Bromocresol Purple LAB-AGAR™ Base	PS 757E	500 g
KF Streptococcal LAB-AGAR™	PS 124	500 g
Kimming Fungi LAB-AGAR™ Base	PS 758E	500 g
King Medium Base	PS 759E	500 g
Kligler Iron LAB-AGAR™	PS 760E	500 g
Kligler Iron LAB-AGAR™ Modified	PS 761E	500 g
Kligler Iron LAB-AGAR™	PS 188	500 g
Kligler Iron LAB-AGAR™ acc.to ISO 10273	PS 08	500 g
Kohn Medium No. 1	PS 762E	500 g
Kohn Medium No.2	PS 763E	500 g
Koser Citrate Broth	PS 30	500 g
Koser Citrate Medium	PS 764E	500 g
L.S. Differential Medium Base	PS 765E	500 g
Lactic Acid Bacteria LAB-AGAR™ Base	PS 766E	500 g

Dehydrated culture media
Millieux de culture deshydrates

Description	Cat.No	Package
Lactic Acid Bacteria Selective Broth Base	PS 767E	500 g
Lactose Broth	PS 35	500 g
Lactose Broth	PS 768E	500 g
Lactose Fluid Medium acc. to USP, EP,BP	PS 770E	500 g
Lactose Gelatin Medium Modified	PS 769E	500 g
Lactose Sulfite Broth Base	PS 162	500 g
Lauryl Sulfate Broth	PS 56	500 g
Lauryl Sulphate LAB-AGAR™	PS 116	500 g
Lauryl Sulphate Tryptose Broth Base, Modified	PS 772E	500 g
Lauryl Tryptose Broth	PS 771E	500 g
Lauryl Tryptose Mannitol Tryptophan Broth	PS 773E	500 g
LB Broth	PS 31	500 g
LB LAB-AGAR™	PS 09	500 g
LEB Broth	PS 420	500 g
LEE Broth	PS 283	500 g
Lees Multidifferential LAB-AGAR™	PS 774E	500 g
Legionella B.C.Y.E. LAB-AGAR™ Base	PS 775E	500 g
Legionella CYE LAB-AGAR™ Base	PS 163	500 g
Lethen Broth	PS 777E	500 g
Lethen Broth Modified	PS 778E	500 g
Lethen LAB-AGAR™ Modified	PS 776E	500 g
Lethen Modified LAB-AGAR™ Base	PS 152	500 g
Lethen w. Triton X100 Broth	PS 779E	500 g
Leuconostoc LAB-AGAR™	PS 247	500 g
Levine`a LAB-AGAR™	PS 235	500 g
Lima Bean LAB-AGAR™	PS 1049E	500 g
Listeria acc.to Oxford LAB-AGAR™ Base	PS 104	500 g
Listeria acc.to Palcam LAB-AGAR™ Base	PS 102	500 g
Listeria Enrichment acc. to FDA Broth	PS 781E	500 g
Listeria Enrichment acc.to Lovett Broth Base (TSYEB)	PS 209	500 g
Listeria Enrichment UVM Broth Base	PS 210	500 g
Listeria Enrichment UVM Medium Base	PS 782E	500 g
Listeria Motility Medium	PS 783E	500 g
Listeria Oxford Medium Base	PS 784E	500 g
Listeria Palcam Lab-AGAR™ Base	PS 780E	500 g
Listeria Selective acc. FDA Broth Base	PS 786E	500 g
Listeria Selective Twin Pack LAB-AGAR™	PS 785E	500 g
Litmus Milk Medium	PS 787E	500 g
Liver Broth	PS 788E	500 g
Liver Broth	PS 129	500 g
Liver Meat LAB-AGAR™	PS 789E	500 g
Lowenstein Jensen Medium Base	PS 790E	500 g
Lowenstein-Jensen Medium Base	PS 241	500 g
LPB Medium	PS 588	500 g
LPT Dilution Broth	PS 302	500 g
LT100 Eugon Broth	PS 310	500 g
Luria Bertani Broth, Miller	PS 792E	500 g
Luria Bertani LAB-AGAR™ Miller	PS 791E	500 g
Lyngby Iron LAB-AGAR™	PS 273	500 g
Lysine Decarboxylase Broth	PS 793E	500 g

Dehydrated culture media

Millieux de culture deshydrates

Description	Cat.No	Package
Lysine Decarboxylase Broth	PS 150	500 g
Lysine Decarboxylase Broth acc.to ISO 6579	PS 57	500 g
Lysine Iron LAB-AGAR™	PS 794E	500 g
Lysine Iron LAB-AGAR™	PS 167	500 g
Lysine Medium Base	PS 795E	500 g
M 17 acc. to Terzaghi Broth	PS 798E	500 g
M 17 acc. to Terzaghi LAB-AGAR™	PS 796E	500 g
M 17 w/o Lactose LAB-AGAR™	PS 797E	500 g
M-17 acc. to Terzaghi LAB-AGAR™	PS 103	500 g
M-17 Broth	PS 108	500 g
Mac Conkey Broth	PS 12	500 g
Mac Conkey LAB-AGAR™	PS 10	500 g
Mac Conkey Nr 1 LAB-AGAR™ acc.to USP	PS 183	500 g
Mac Conkey Nr 2 LAB-AGAR™	PS 11	500 g
Mac Conkey Purple w. Taurocholate Broth	PS 802E	500 g
Mac Conkey Purple w/o Taurocholate Broth	PS 803E	500 g
Mac Conkey Sorbitol LAB-AGAR™ Base	PS 805E	500 g
Mac Conkey w. Neutral Red Broth	PS 804E	500 g
Mac Conkey w/o crystal violet and w/o sodium chloride LAB-AGAR™	PS 111	500 g
Mac Conkey w/o CV LAB-AGAR™	PS 800E	500 g
Mac Conkey w/o CV LAB-AGAR™ acc. USP	PS 801E	500 g
Malachite Green Broth	PS 806E	500 g
Malonate Broth	PS 807E	500 g
Malonate Broth Ewing Modified	PS 808E	500 g
Malt Extract Broth	PS 78	500 g
Malt Extract Broth	PS 810E	500 g
Malt Extract LAB-AGAR™ Base	PS 809E	500 g
Malt Extract LAB-AGAR™	PS 77	500 g
Malt Extract Peptone LAB-AGAR™	PS 260	500 g
Mannitol Lysine LAB-AGAR™	PS 811E	500 g
Mannitol Salt acc.to Chapman LAB-AGAR™ (acc.to USP)	PS 186	500 g
Mannitol Salt Broth	PS 256	500 g
Mannitol Salt LAB-AGAR™	PS 13	500 g
Mannitol Selenite Broth	PS 813E	500 g
Maximum Recovery Diluent	PS 814E	500 g
M-Broth	PS 815E	500 g
Mc Bridge Listeria LAB-AGAR™ Base	PS 816E	500 g
M-CP LAB-AGAR™ Base	PS 817E	500 g
m-CP LAB-AGAR™ Base	PS 138	500 g
MEA LAB-AGAR™ BASE acc. to ISO 16212	PS 232	500 g
M-Endo Broth	PS 819E	500 g
M-Endo LES LAB-AGAR™	PS 818E	500 g
M-Enterococcus LAB-AGAR™ Base	PS 820E	500 g
M-FC Broth Base	PS 822E	500 g
M-FC LAB-AGAR™ Base	PS 821E	500 g
MICROBIAL CONTENT TEST LAB-AGAR™	PS 238	500 g
Middlebrook 7H9 Broth Base	PS 824E	500 g
Middlebrook 7H9 LAB-AGAR™ Base	PS 823E	500 g
MIL Medium	PS 825E	500 g
Milk LAB-AGAR™	PS 826E	500 g

Dehydrated culture media
Millieux de culture deshydrates

Description	Cat.No	Package
Milk w. Cetrimide Twin Pack LAB-AGAR™	PS 827E	500 g
Mineral Glutamate Modified Medium	PS 323	500 g
Minerals Glutamate Twin Pack LAB-AGAR™ Base	PS 828E	500 g
MIO LAB-AGAR™	PS 58	500 g
Mitis salivarius LAB-AGAR™ Base	PS 830E	500 g
MIU Medium Base	PS 831E	500 g
M-Lauryl Sulphate Broth	PS 832E	500 g
m-LST Broth Base	PS198	500 g
m-LST Broth Base	PS 198	500 g
M-MacConkey Broth	PS 833E	500 g
Mossel EE Broth	PS 86	500 g
Motility Indole Ornithine Medium	PS 829E	500 g
Motility LAB-AGAR™ acc.to ISO 11290	PS 536	500 g
Motility LAB-AGAR™ acc.to ISO 6579	PS 516	500 g
Motility Nitrate Medium Base	PS 115	500 g
Motility Nitrate Medium Buffered	PS 834E	500 g
Motility Test Medium	PS 835E	500 g
MRS Broth	PS 60	500 g
MRS LAB-AGAR pH 5,7	PS 271	500 g
MRS LAB-AGAR™	PS 59	500 g
MRS Lactobacillus Broth	PS 838E	500 g
MRS Lactobacillus LAB-AGAR™	PS 836E	500 g
MRS Lactobacillus LAB-AGAR™ Modified	PS 837E	500 g
MRVP Medium	PS 92	500 g
MR-VP Medium acc. to ISO 6579	PS 839E	500 g
MSRV Medium	PS 225	500 g
MSRV Medium Base	PS 203	500 g
m-TSB Broth (w. novobiocin)	PS 550	500 g
m-TSB Broth Base	PS 197	500 g
MUD SF Broth Base	PS 840E	500 g
Mueller Hinton 2 LAB-AGAR™	PS 79	500 g
Mueller Hinton Broth	PS 15	500 g
Mueller Hinton Broth	PS 842E	500 g
Mueller Kauffman Broth Base	PS 84	500 g
Mueller Kaufman w. Novobiocin and Brilliant Green Broth Base acc. to ISO 6579	PS 147	500 g
MUG EC Broth Modified	PS 844E	500 g
MUG Lauryl Sulphate Broth Modified	PS 845E	500 g
MUG Mac Conkey LAB-AGAR™	PS 846E	500 g
MUG Plate Count LAB-AGAR™	PS 847E	500 g
MUG Violet Red Bile LAB-AGAR™	PS 848E	500 g
Mycobiotic LAB-AGAR™	PS 161	500 g
Mycoplasma PPLO Broth Base	PS 850E	500 g
Mycoplasma PPLO LAB-AGAR™ Base	PS 849E	500 g
MYP acc. to Mossel LAB-AGAR™ Base Modified	PS 851E	500 g
NBB Broth Base	PS 365	500 g
NBB LAB-AGAR™ Base	PS 410	500 g
Neutralizing Fluid Medium	PS 852E	500 g
Niacin Assay Medium	PS 853E	500 g
Nitrate Broth	PS 854E	500 g
Nitrate LAB-AGAR™	PS 1040E	500 g

Dehydrated culture media

Milieux de culture deshydrates

Description	Cat.No	Package
NOCIVE BREWERS BACTERIA AGAR BASE MODIFIED	PS 315	500 g
NOCIVE BREWERS BACTERIA BROTH BASE MODIFIED	PS 480	500 g
Nocive Brewers Bacteria Broth Base Modified	PS 480	500 g
Nocive Brewers Bacteria LAB-AGAR™ Base Modified	PS 315	500 g
Nutrient 1% Peptone Broth	PS 865E	500 g
Nutrient at pH 6,8 LAB-AGAR™	PS 859E	500 g
Nutrient at pH 6,9 Broth	PS 864E	500 g
Nutrient at pH 7,0 LAB-GAR™	PS 860E	500 g
Nutrient Broth	PS 90	500 g
Nutrient Broth	PS 863E	500 g
Nutrient Broth No 2	PS 204	500 g
Nutrient E LAB-AGAR™	PS 330	500 g
Nutrient Gelatin	PS 61	500 g
Nutrient Gelatin	PS 866E	500 g
Nutrient Glucose Broth	PS 507	500 g
Nutrient LAB-AGAR™	PS 855E	500 g
Nutrient LAB-AGAR™ acc. ISO 13270	PS 856E	500 g
Nutrient LAB-AGAR™	PS 85	500 g
Nutrient LAB-AGAR™ acc.to ISO 21528	PS 166	500 g
Nutrient LAB-AGAR™ acc.to ISO 9308-1	PS 194	500 g
Nutrient LAB-AGAR™ acc.to EN 12780; ISO 16266	PS 156	500 g
Nutrient No 2 LAB-AGAR™	PS 858E	500 g
Nutrient Oxidase LAB-AGAR™	PS 857E	500 g
Nutrient S LAB-AGAR™	PS 307	500 g
Nutrient Semisolid LAB-AGAR™	PS 861E	500 g
Nutrient w. Manganese LAB-AGAR™	PS 862E	500 g
O.G.A. LAB-AGAR™	PS 126	500 g
OF Basal Medium (Hugh Leifson)	PS 75	500 g
OF Basal Staphylococcus Medium	PS 867E	500 g
Orange Serum LAB-AGAR™	PS 868E	500 g
Orange Serum LAB-AGAR™	PS 265	500 g
Ornithine Decarboxylase Broth	PS 869E	500 g
ORSIM LAB-AGAR™ Base	PS 286	500 g
Osmophilic LAB-AGAR™	PS 157	500 g
Oxacillin Resistance Screening LAB-AGAR™ Base	PS 870E	500 g
Oxytetracycline LAB-AGAR™ Base	PS 277	500 g
Oxytetracycline Glucose Yeast Extract LAB-AGAR™ Base	PS 871E	500 g
Page'a Buffer Solution	PS 299	500 g
Pantothenate Assay Medium	PS 882E	500 g
Park and Sanders Broth Base	PS 883E	500 g
PCA LAB-AGAR™ acc. to ISO 4833-2	PS 261	500 g
PE-2 Medium	PS 1041E	500 g
PEMBA LAB-AGAR™ Base	PS 1032E	500 g
PEMBA Selective LAB-AGAR™ Base	PS 180	500 g
Peptone Water	PS 127	500 g
Peptone Water	PS 885E	500 g
Peptone Water w. Andrade Indicator	PS 887E	500 g
Peptone Water w. Phenol Red	PS 886E	500 g
Peptone Water with Lactose	PS 105	500 g
Perfringens (OPSP) LAB-AGAR™ Base	PS 889E	500 g

Dehydrated culture media
Millieux de culture deshydrates

Description	Cat.No	Package
Perfringens LAB-AGAR TM Base	PS 888E	500 g
pH 7,0 Pharmacopeia Diluent	PS 181	500 g
Phenol Red Broth Base	PS 890E	500 g
Phenol Red Broth Base	PS 191	500 g
Phenol Red Lactose Broth	PS 891E	500 g
Phenol Red Sucrose Broth	PS 1035E	500 g
Phenylalanine LAB-AGAR TM	PS 892E	500 g
Phenylalanine LAB-AGAR TM	PS 62	500 g
Phenylalanine Malonate Broth	PS 893E	500 g
Phenylketonuria (PKU) Test LAB-AGAR TM Base	PS 895E	500 g
Phosphate Buffer acc. to ISO 16654	PS 894E	500 g
Pirazynamidase LAB-AGAR TM	PS 538	500 g
Plate Count (Standard Methods) LAB-AGAR TM	PS 896E	500 g
Plate Count Tryptone Glucose Yeast Extract LAB-AGAR TM	PS 897E	500 g
PLET B. anthracis LAB-AGAR TM Base	PS 898E	500 g
Potassium Cyanide w/o KCN Broth Base	PS 899E	500 g
Potato Carrot LAB-AGAR TM	PS 1044E	500 g
Potato Dextrose Broth	PS 231	500 g
Potato Dextrose LAB-AGAR TM acc. to USP	PS 900E	500 g
Potato Dextrose LAB-AGAR TM	PS 107	500 g
Potato Dextrose Sucrose LAB-AGAR TM	PS 1048E	500 g
PPLO w/o cristal violet LAB-AGAR TM Base	PS 140	500 g
PPLO w/o cristal violet Broth Base	PS 141	500 g
Preston Broth Base	PS 255	500 g
Preston LAB-AGAR TM Base	PS 901E	500 g
PSB Broth Modified	PS 902E	500 g
Pseudomonas CFC LAB-AGAR TM Base	PS 903E	500 g
Pseudomonas CN LAB-AGAR TM Base	PS 158	500 g
Pseudomonas F LAB-AGAR TM (KING B)	PS 64	500 g
Pseudomonas Fluorescein LAB-AGAR TM Base	PS 904E	500 g
Pseudomonas Isolation LAB-AGAR TM (acc. to USP)	PS 184	500 g
Pseudomonas LAB-AGAR TM acc. to ISO EN 12780, ISO 16266	PS 189	500 g
Pseudomonas LAB-AGAR TM Base (CFC)	PS 144	500 g
Pseudomonas P LAB-AGAR TM - King A	PS 592	500 g
Pseudomonas Pyocyanin LAB-AGAR TM Base	PS 905E	500 g
Purple Broth Base acc. ISO 11290	PS 906E	500 g
R2 LAB-AGAR TM (acc.to USP)	PS 185	500 g
R2A Broth	PS 1042E	500 g
R2A LAB-AGAR TM	PS 130	500 g
Raka-Ray LAB-AGAR TM Base	PS 120	500 g
Rappaport Soy Broth	PS 65	500 g
Rappaport Vassiliadis Broh acc. to European Pharmacopoeia	PS 250	500 g
Rappaport Vassiliadis Broth acc. ISO 6340	PS 919E	500 g
Rappaport Vassiliadis Medium	PS 907E	500 g
Rappaport Vassiliadis Medium Modified	PS 908E	500 g
Rappaport Vassiliadis RVSM Broth	PS 909E	500 g
Rappaport Vassiliadis Soya Broth	PS 920E	500 g
Reincorfed Clostridial (RCA) LAB-AGAR TM	PS 910E	500 g
Reincorfed Clostridial LAB-AGAR TM	PS 142	500 g
Reincorfed Clostridial Medium	PS 143	500 g

Dehydrated culture media

Millieux de culture deshydrates

Description	Cat.No	Package
Renicorfed Clostridial (RCM)Broth	PS 911E	500 g
Rice Extract LAB-AGAR™	PS 912E	500 g
Ringer Salt Solution Powder	PS 913E	500 g
Rogosa Broth	PS 915E	500 g
Rogosa LAB-AGAR™	PS 914E	500 g
Rogosa SL LAB-AGAR™	PS 916E	500 g
Rose Bengal Chloramphenicol LAB-AGAR™	PS 918E	500 g
Rose Bengal LAB-AGAR™ Base	PS 917E	500 g
Rose Bengal LAB-AGAR™	PS 66	500 g
RPMI MOPS LAB-AGAR™ BASE	PS 233	500 g
S. S. LAB-AGAR™	PS 921E	500 g
S. S. LAB-AGAR™ Modified	PS 922E	500 g
Sabouraud 2% Dex. w. Chloramphenicol LAB-AGAR™	PS 227	500 g
Sabouraud 2% Dextrose LAB-AGAR™ Base Modified	PS 923E	500 g
Sabouraud 2% Maltose Broth	PS 924E	500 g
Sabouraud 4% Maltose LAB-AGAR™	PS 925E	500 g
Sabouraud Dextrose 2% LAB-AGAR™	PS 547	500 g
Sabouraud Dextrose Broth	PS 146	500 g
Sabouraud Dextrose LAB-AGAR™	PS 16	500 g
Sabouraud Dextrose w. Chloramphenicol LAB-AGAR™	PS 32	500 g
Sabouraud Dextrose w. Chloramphenicol LAB-AGAR™ (acc. to EU. PHARM)	PS 192	500 g
Sabouraud Dextrose w.Chloramphenicol and Cycloheximide LAB-AGAR™	PS 94	500 g
Sabouraud Dextrose w.Chloramphenicol and Gentamycine LAB-AGAR™	PS 248	500 g
Sabouraud Fluid Medium	PS 311	500 g
Sabouraud Fluid Medium	PS 930E	500 g
Sabouraud Maltose Broth	PS 229	500 g
Sabouraud Maltose LAB-AGAR™	PS 228	500 g
Saline Buffered (PBS) Solution	PS 269	500 g
Saline Peptone Water	PS 118	500 g
Salmonella / Shigella SSDC LAB-AGAR™	PS 170	500 g
Salmonella acc. to ÖNÖZ LAB-AGAR™	PS 931E	500 g
SALMONELLA DIFFERENTIAL LAB-AGAR™ BASE (Rambach)	PS 240	500 g
Salmonella Differential Twin Pack LAB-AGAR™ Modified	PS 932E	500 g
Salmonella Shigella LAB-AGAR™	PS 17	500 g
Salt Meat Broth	PS 933E	500 g
Salt Polymixin Broth Base	PS 934E	500 g
Schaedler Agar / Schaedler LAB-AGAR™	PS 935E	500 g
Schaedler Broth	PS 19	500 g
Schaedler Broth	PS 936E	500 g
Schaedler LAB-AGAR™	PS 18	500 g
SDA/LMDA	PS 317	500 g
Selenite Broth Base	PS 937E	500 g
Selenite Cystine Broth	PS 36	500 g
Selenite Cystine Twin Pack Broth	PS 938E	500 g
Selenite Cystine Twin Pack Broth	PS 942E	500 g
Selenite Cystine Twin Pack Medium acc.to USP	PS 939E	500 g
Selenite F Twin Pack Broth	PS 941E	500 g
Semisolid LAB-AGAR™ acc.to ISO 6579	PS 943E	500 g
Semisolid Rappaport Vassiliadis Medium, Modified	PS 944E	500 g
Sensitivity Test LAB-AGAR™	PS 741E	500 g

Dehydrated culture media
Millieux de culture deshydrates

Description	Cat.No	Package
Sensitivity Test No. 2 LAB-AGAR™	PS 742E	500 g
SFP LAB-AGAR™ Base	PS 945E	500 g
Shapton Medium	PS 946E	500 g
SIM Medium	PS 947E	500 g
Simmons Citrate LAB-AGAR™	PS 948E	500 g
Simmons Citrate LAB-AGAR™	PS 20	500 g
Skim Milk LAB-AGAR™ Modified	PS 949E	500 g
Skim Milk Plate Count LAB-AGAR™	PS 950E	500 g
Slanetz and Bartley LAB-AGAR™	PS 67	500 g
Slanetz and Bartley Medium	PS 951E	500 g
Sodium Selenite Broth	PS 68	500 g
Sorbitol Mac Conkey LAB-AGAR™	PS 21	500 g
Sorbitol Peptone with Bile Salts (PSB)	PS 63	500 g
Soy Tryptone LAB-AGAR™	PS 133	500 g
Soyabean Bile Broth Base Modified	PS 952E	500 g
Soyabean Bile w. Novobiocin Broth acc. ISO 16654	PS 953E	500 g
SPS acc.to Angelott LAB-AGAR™	PS 954E	500 g
Standard Methods (Plate Count Agar) LAB-AGAR™	PS 37	500 g
Standard Methods w. powdered Milk LAB-AGAR™	PS 38	500 g
Standard Nutrient No 1 Broth	PS 956E	500 g
Standard Nutrient No 1 LAB-AGAR™	PS 955E	500 g
Staphylococcus No 110 LAB-AGAR™	PS 957E	500 g
Starch LAB-AGAR™ Agar	PS 1033E	500 g
Streptococcus agalactiae LAB-AGAR™ Base	PS 958E	500 g
Streptococcus Enrichment (SE) Broth	PS 959E	500 g
Streptococcus Selective Broth	PS 70	500 g
Streptococcus Selective LAB-AGAR™	PS 69	500 g
Stuart Transport Medium	PS 981E	500 g
Sulphate API w/o Sodium Lactate LAB-AGAR™	PS 960E	500 g
Synthetic Sea Salt	PS 961E	500 g
T.A.T.Broth Base	PS 962E	500 g
TBX Tryptone Bile X-glucuronide LAB-AGAR™	PS 135	500 g
TCBS Cholera Medium	PS 964E	500 g
TCBS LAB-AGAR™	PS 963E	500 g
TCBS LAB-AGAR™	PS 71	500 g
Teepol Twin Pack Broth Modified	PS 965E	500 g
Tergitol 7 LAB-AGAR™ Base Modified	PS 966E	500 g
Tetrathionate Brilliant Green Bile Broth	PS 967E	500 g
Tetrathionate Broth Base acc.to Mueller Kauffman	PS 969E	500 g
Tetrathionate Crystal Violet Enrichment Broth	PS 971E	500 g
Tetrathionate Mueller Kauffman Broth Base acc. to ISO 3565	PS 970E	500 g
Tetrathionate w/o Brilliant Green Broth Base	PS 968E	500 g
Thayer Martin Medium Base	PS 972E	500 g
Thermoacidurans LAB-AGAR™	PS 1043E	500 g
Thioglycollate Broth acc. to NIH	PS 974E	500 g
Thioglycollate Fluid Medium	PS 43	500 g
Thioglycollate Fluid Medium acc.to USP	PS 973E	500 g
Thioglycollate Fluid w. Resazurin Medium	PS 976E	500 g
Thioglycollate w/o Indicator Medium	PS 977E	500 g
Thioglycollate Enriched Broth	PS 975E	500 g

Dehydrated culture media

Milieux de culture deshydrates

Description	Cat.No	Package
Tinsdale LAB-AGAR™ Base	PS 278	500 g
Tinsdale Agar Base/ Tinsdale LAB-AGAR™ Base	PS 978E	500 g
Todd Hewitt Broth	PS 72	500 g
Tomato Juice LAB-AGAR™	PS 980E	500 g
Tood Hewitt Broth	PS 979E	500 g
TOS MEDIUM Base	PS 300	500 g
TPY Broth	PS 534	500 g
Tributylin LAB-AGAR™ Base	PS 982E	500 g
Trichomonas LAB-AGAR™ Base	PS 982E	500 g
Trichomonas Medium Base	PS 175	500 g
Triple Sugar Iron LAB-AGAR™	PS 984E	500 g
Triple Sugar Iron LAB-AGAR™	PS 44	500 g
Trypticasein Soy Broth (TSB)	PS 23	500 g
Trypticasein Soy LAB-AGAR™ (TSA)	PS 22	500 g
Tryptone Bile LAB-AGAR™	PS 985E	500 g
Tryptone Glucose Beef Extract (TGBEA) LAB-AGAR™	PS 987E	500 g
Tryptone Glucose Extract (TGYES) LAB-AGAR™	PS 988E	500 g
Tryptone Glucose Extract LAB-AGAR™	PS 505	500 g
Tryptone Purple Yeast Extract LAB-AGAR™	PS 998E	500 g
Tryptone Salt Broth	PS 989E	500 g
Tryptone Soya Broth	PS 992E	500 g
Tryptone Soya + Lecithin + Tween 80 LAB-AGAR™	PS 991E	500 g
Tryptone Soya LAB-AGAR™	PS 990E	500 g
Tryptone Soya Yeast Extract (TSYEA) LAB-AGAR™	PS 993E	500 g
Tryptone Soya Yeast Extract (TSYEB) Broth	PS 994E	500 g
Tryptone Sucrose Tetrazolium (TSTA) LAB-AGAR™ Base	PS 995E	500 g
Tryptone Water	PS 201	500 g
Tryptone Water acc. to ISO 7251	PS 996E	500 g
Tryptone with Bile Salts LAB-AGAR™ (TBA)	PS 128	500 g
Tryptone Yeast Extract Cystine Sucrose LAB-AGAR™ Base	PS 999E	500 g
Tryptone Yeast Extract LAB-AGAR™	PS 997E	500 g
Tryptophan Culture Broth	PS 132	500 g
Tryptophan Medium	PS 1000E	500 g
Tryptose Blood LAB-AGAR™ Base	PS 1002E	500 g
Tryptose Broth	PS 1003E	500 g
Tryptose LAB-AGAR™	PS 1001E	500 g
Tryptose Lauryl Sulfate Mannitol Tryptophan Broth	PS 202	500 g
Tryptose Lauryl Sulfate MUG Broth	PS 200	500 g
Tryptose Phosphate Broth	PS 1004E	500 g
Tryptose Sulphite Neomycin LAB-AGAR™	PS 1005E	500 g
TSA + 4N LAB-AGAR™	PS 287	500 g
TSC LAB-AGAR™ Base acc. to ISO 14189	PS 291	500 g
TSC LAB-AGAR™	PS 96	500 g
TSN LAB-AGAR™	PS 91	500 g
TSYEA LAB-AGAR™ acc.to ISO 11290	PS 193	500 g
TTC Tergitol-7 LAB-AGAR™	PS 125	500 g
UBA Universal Beer LAB-AGAR™	PS 264	500 g
Universal Beer Agar (UBA)	PS 293	500 g
Uniwersal Beer (UBA) LAB-AGAR™	PS 1006E	500 g
Urea Broth	PS 1008E	500 g

Dehydrated culture media
Millieux de culture deshydrates

Description	Cat.No	Package
Urea Broth	PS 593	500 g
Urea Broth Base	PS 506	500 g
Urea Broth Base	PS 1009E	500 g
Urea Christensen LAB-AGAR™ Base	PS 1007E	500 g
Urea Indol Broth	PS 134	500 g
UREA LAB-AGAR™ Base	PS 24	500 g
UREA LAB-AGAR™ Base acc.to ISO 6579	PS 211	500 g
V-8 Juice LAB-AGAR™	PS 1045E	500 g
Violet Red Bile Glucose LAB-AGAR™	PS 1010E	500 g
Violet Red Bile Lactose LAB-AGAR™	PS 1011E	500 g
Violet Red Bile with Glucose LAB-AGAR™	PS 73	500 g
Violet Red Bile with Lactose LAB-AGAR™	PS 74	500 g
Vitamin B12 Assay Medium	PS 1012E	500 g
Vogel-Johnson LAB-AGAR	PS 236	500 g
Vogel-Johnson w/o Tellurite LAB-AGAR™ Base	PS 1013E	500 g
V-P Broth Modified	PS 1014E	500 g
VP Medium acc. to ISO 11290	PS 294	500 g
VRBLG LAB-AGAR™	PS 172	500 g
Wilkins Chalgren Anaerobic Broth Base	PS 1016E	500 g
Wilkins Chalgren Anaerobic LAB-AGAR™ Base	PS 1015E	500 g
Wilkins Chalgren LAB-AGAR™	PS 25	500 g
Willis and Hobbs Medium Base	PS 1017E	500 g
Wilson Blair LAB-AGAR™ Base	PS 1018E	500 g
Wilson Blair LAB-AGAR™ Base (dla beztlenowców)	PS 153	500 g
WL Nutrient Broth	PS 1020E	500 g
WL Nutrient LAB-AGAR™	PS 1019E	500 g
WL Nutrient LAB-AGAR™	PS 97	500 g
Wort Broth	PS 218	500 g
Wort Broth	PS 1022E	500 g
Wort LAB-AGAR™	PS 1021E	500 g
Wort LAB-AGAR™	PS 81	500 g
X.L.D. LAB-AGAR™ acc.to ISO 6579	PS 45	500 g
X.L.T.4 LAB-AGAR™ Base	PS 131	500 g
XLT4 LAB-AGAR™ Base	PS 1024E	500 g
Xylose Lysine Desoxycholate LAB-AGAR™ (medium K acc. Pharmacopoeia)	PS 308	500 g
Yeast Extract Dextrose LAB-AGAR™	PS 196	500 g
Yeast Extract LAB-AGAR™	PS 1025E	500 g
Yeast Extract LAB-AGAR™ acc.to ISO 6222	PS 82	500 g
Yeast Malt (YM) LAB-AGAR™	PS 1026E	500 g
Yeast Malt Broth	PS 1027E	500 g
Yeast Mold Broth	PS 1047E	500 g
Yeast Mold LAB-AGAR™	PS 1046E	500 g
Yeast Nitrogen Base	PS 26	500 g
Yersinia Identification Broth Base	PS 1028E	500 g
Yersinia Isolation LAB-AGAR™	PS 1029E	500 g
Yersinia Selective LAB-AGAR™ Base	PS 1030E	500 g
Yersinia Selective LAB-AGAR™ Base	PS 145	500 g
YPG LAB-AGAR™ Base	PS 1031E	500 g

Ingredients of media**Composants des milieux**

Description	Cat.No	Package
AGARS		
BACTERIOLOGICAL LAB -AGAR™	AB 03	500 g
LAB -AGAR- AGAR™	AB 04	500 g
PURIFIED LAB -AGAR™	AB 06	500 g
VITRO LAB -AGAR™	AB 08	500 g
PLANT PROPAGATION LAB -AGAR™-	AB 12	500 g
PLANT TC LAB-AGAR™	AB 14	500 g
PHARMACEUTICAL LAB -AGAR™	AB 16	500 g
PEPTONES		
Bacteriological Peptone	PB 01	500 g
Tryptone	PB 02	500 g
Beef Extract	PB 03	500 g
Casein Peptone	PB 04	500 g
Casein Peptone H	PB 05	500 g
Yeast Extract Mikrogranulat	PB 06	500 g
Dehydrated Oxbile	PB 07	500 g
Gelatin Peptone	PB 08	500 g
Meat Peptone	PB 09	500 g
Polypeptone	PB 10	500 g
Proteose Peptone	PB 11	500 g
Bile Salts Nr 3	PB 12	500 g
Soy Peptone	PB 13	500 g
Tryptose	PB 14	500 g
Heart Infusion	PB 15	500 g
Lactoalbumin hydrolysate	PB 16	500 g
Malt Extract	PB 17	500 g
Proteose Peptone Nr 3	PB 18	500 g
Mycological Peptone	PB 40	500 g
Gelatin Powdered	PB 42	500 g

Sterile animal blood in bottles**Sang sterile animalier en bouteille**

Description	Cat.No	Package
Sterile defibrinated Sheep Blood	SL 0160	1 x 100 ml
Sterile defibrinated Sheep Blood	SL 0160-200	1 x 200 ml
Sterile defibrinated Sheep Blood	SL 0160-500	1 x 500 ml
Sterile defibrinated Horse Blood	SL 0150	1 x 100 ml
Sterile defibrinated Horse Blood	SL 0150-200	1 x 200 ml
Sterile defibrinated Horse Blood	SL 0150-500	1 x 500 ml

Supplements

Supplements

Description	Cat.No	Package size
Ampicillin Supplement	SL 0100	10 vials 1/500 ml
Bacillus cereus Supplement	SL 0002	10 vials 1/500 ml
Blaser-Wang Supplement	SL 0095	10 vials 1/500 ml
Bolton Selective Supplement	SL 0050	10 vials 1/500 ml
Bordetella Supplement	SL 0003	10 vials 1/500 ml
Brucella Supplement	SL 0004	10 vials 1/500 ml
Burkholderia cepacia Supplement	SL 0068	10 vials 1/500 ml
Butzler Modified Supplement	SL 0083	10 vials 1/500 ml
Butzler Supplement	SL 0089	10 vials 1/500 ml
Campylobacter CCDA Selective Supplement	SL 0005	10 vials 1/500 ml
Campylobacter SKIRROW Supplement	SL 0006	10 vials 1/500 ml
Cefoxitin Supplement	SL 0063	10 vials 1/500 ml
Cefsulodin Novobiocin Supplement	SL 0066	10 vials 1/500 ml
Cefsulodin Vancomycin Supplement	SL 0067	10 vials 1/500 ml
CFC Pseudomonas Supplement	SL 0007	10 vials 1/500 ml
Chloramphenicol Supplement	SL 0058	10 vials 1/500 ml
Chromogenic Listeria Kit Supplement acc.to ISO 11290	SL 0001	2x5 vials 1+1/500 ml
Clostridium difficile Selective Supplement	SL 0048	10 vials 1/500 ml
CN Supplement	SL 0057	10 vials 1/500 ml
CNA Supplement	SL 0009	10 vials 1/500 ml
COA Supplement	SL 0065	10 vials 1/500 ml
C-T Supplement	SL 0056	10 vials 1/500 ml
Dermatophytes Supplement	SL 0064	10 vials 1/500 ml
Dermatophytes Modified Supplement	SL 0080	10 vials 1/500 ml
Egg's Yolk Polimixin Sterile Emulsion	SL 0037	100 ml
Egg's Yolk Sterile Emulsion	SL 0035	100 ml
Egg's Yolk Tellurite Sterile Emulsion	SL 0036	100 ml
Endo Basic Fuchsin Solution 10%	SL 0010	60 ml
Ferric Ammonium Citrate Supplement	SL 0011	10 vials 1/500 ml
Ferric Ammonium Citrate Supplement	SL 0011B	10 vials 1/225ml
Ferric Ammonium Citrate Supplement	SL 0062	10 vials 1/2,25L
Fraser Listeria Selective Supplement	SL 0012	10 vials 1/500 ml
Gardnerella Selective Supplement	SL 0082	10 vials 1/500 ml
Glycerol	GLC 80500	500 ml
Haemophilus Test Supplement	SL 0054	10 vials 1/500 ml
Half Fraser Listeria Selective Supplement	SL 0014	10 vials 1/500 ml
Half Fraser Listeria Selective Supplement	SL 0108	10 vials 1/2,25 L
Half Fraser Listeria Supplement	SL 0107	10 vials 1/225ml
Horse Serum	SL 0116	100 ml
ITC Selective Supplement	SL 0015	10 vials 1/500 ml
Karmali Modified Selective Supplement	SL 0118	10 vials 1/500 ml
Karmali Selective Supplement	SL 0049	10 vials 1/500 ml
L.C.A.T. Supplement	SL 0016	10 vials 1/500 ml
10% LACTIC ACID	SL 0079	10 vials x1 ml
10% Lactic Acid	SL 0086	100 ml
Lecithine	SL 0500	500 g
Legionella BCYE without cysteine Supplement	SL 0046	10 vials 1/100 ml
Legionella BCYE without cysteine Supplement	SL 0047	10 vials 1/500 ml
Legionella CYE Growth Supplement	SL 0017	10 vials 1/100 ml

Supplements

Supplements

Description	Cat.No	Package size
Legionella CYE Growth Supplement	SL 0053	10 vials 1/500 ml
Legionella GVPC Supplement	SL 0018	10 vials 1/500 ml
Legionella Selective AB Supplement	SL 0087	10 vials 1/100 ml
Legionella Selective AB Supplement	SL 0096	10 vials 1/500 ml
LEGIONELLA SELECTIVE MWY SUPPLEMENT	SL 0020	10 vials 1/100 ml
LEGIONELLA SELECTIVE MWY SUPPLEMENT	SL 0085	10 vials 1/500 ml
Listeria Modified acc. to Oxford Supplement	SL 0081	10 vials 1/500 ml
Listeria Selective acc. to Oxford Supplement	SL 0021	10 vials 1/500 ml
Listeria Selective acc. to Palcam Supplement	SL 0022	10 vials 1/500 ml
Listeria Selective UVM Supplement	SL 0024	10 vials 1/500 ml
Listeria Selective UVM II Supplement	SL 0025	10 vials 1/500 ml
L-Malic Acid Supplement	SL 0088	10 vials 1/500 ml
Lysed Horse Blood	SL 0103	100 ml
m-CP Selective Supplement	SL 0026	2x5 vials 1+1/500 ml
MSRV Selective Supplement	SL 0027	10 vials 1/1L
MSRV Selective Supplement	SL 0027-5	10 vials 1/500 ml
MUP Selective Supplement/ TOS Supplement	SL0076	10 vials 1/500 ml
NAD Supplement	SL 0307	10 vials 1/500 ml
Nessler Reagent	NSR 60010	10 ml
Nowobiocin Supplement	SL 0055	10 vials 1/500 ml
O.G.A. Supplement	SL 0028	10 vials 1/500 ml
ORSIM Supplement	SL 0077	10 vials 1/500 ml
Polienrichment Supplement (eq. PolyviteX)	SL 0029	2x5 vials 1+1/500 ml
Potassium Tellurite Sterile Solution 1%	SL 0031/100	100 ml
Potassium Tellurite Sterile Solution 1%	SL 0031/30	30 ml
Potassium Tellurite Sterile Solution 3,5%	SL 0030/100	100 ml
Potassium Tellurite Sterile Solution 3,5%	SL 0030 / 30	30 ml
Potassium Tellurite Sterile Solution 3,5%	SL 0102	10 x 2ml
Preston Campylobacter Supplement	SL 0032	10 vials 1/500 ml
Preston Modified Supplement	SL 0084	10 vials 1/500 ml
Ringer Tablets	SL 0112	100 tablets
RPF Supplement	SL 0034	10 vials 1/100 ml
RPMI-MOPS Supplement	SL 0060	100 ml
Salmonella Chromogenic Supplement 5+5	SL 0061	2x5 vials 1+1/500 ml
Skim Milk Powder	SL 0113	500 g
Sterile Paraffin	SL 0138	100 ml
Sulphamandelate Selective Supplement	SL 0038	10 vials 1/500 ml
Artificial sea water acc. ISO 9308-3	SL 0700-1	100 g
Artificial sea water acc. ISO 9308-3	SL 0700-5	500 g
Tinsdale Enrichment Supplement	SL 0074	10 vials 1/100 ml
Trichomonas Selective Supplement	SL 0039	10 vials 1/500 ml
TRITON X-100 Supplement	SL 0091	100 ml
TSC Supplement	SL 0008	10 vials 1/500 ml
TTC Sterile Solution 1%	SL 0040/100	100 ml
TTC Sterile Solution 5%	SL 0115	100 ml
Tween 80 / 100 ml	TWS 80100	100 ml
Tween 80 / 500 ml	TWS 80500	500 ml
Urea Sterile Solution 40%	SL 0041	100 ml
Urea Sterile Solution 40%	SL 0111	10 x 5ml

Supplements

Supplements

Description	Cat.No	Package size
V.C.A.T. Supplement	SL 0042	10 vials 1/500 ml
V.C.N. Supplement	SL 0043	10 vials 1/500 ml
Vancomycin Supplement	SL 0059	10 vials 1/500 ml
VCC Supplement	SL 0104	10 vials 1/2.25L
XLT 4 Supplement	SL 0045	46ml
Yersinia CIN Selective Supplement	SL 0044	10 vials 1/500 ml

Antibiotic discs

Disques antibiogramme

The Antibiotic Paper Discs are manufactured in a modern facility equipped with automated production lines based on the latest technologies. Each delivery is accompanied by a certificate of quality control. The discs are packed in airtight packagings, type of tube or blister with a desiccant. Production is carried out in accordance with Quality Management System ISO 9001 and ISO 13485. The products are registered and CE marked.

Pour la fabrication des disques antibiogrammes nous utilisons une installation moderne équipée des lignes de production automatiques de dernière génération. Chaque livraison est accompagnée d'un certificat de qualité. Les disques sont emballés en paquets étanches à l'air. Les paquets sont de type tube ou blister ayant un agent déshydratant. La production est conforme au Système de Gestion de la Qualité ISO 9001 et ISO 13485. Les produits sont enregistrés et marqués CE.

Antibiotic paper discs

Disques antibiogramme

Description	Cat.No	Package size
Amikacin 10	E110123	1 x 50 discs
Amikacin 30	E110124	1 x 50 discs
Amoxycillin 10	E110135	1 x 50 discs
Amoxycillin 25	E110132	1 x 50 discs
Amoxycillin/Clavulan acid 2/1	E110117	1 x 50 discs
Amoxycillin/Clavulan acid 20/10	E110116	1 x 50 discs
Ampicillin 10	E110122	1 x 50 discs
Ampicillin 2	E110129	1 x 50 discs
Ampicillin 25	E110125	1 x 50 discs
Ampicillin/Sulbactam 10/10	E110127	1 x 50 discs
Ampicillin/Sulbactam 20/10	E110128	1 x 50 discs
Azithromycin 15	E110145	1 x 50 discs
Aztreonam 30	E110146	1 x 50 discs
Bacitracin 10	E110312	1 x 50 discs
Benzylpenicillin 1	E111606	1 x 50 discs
Carbenicillin 100	E113017	1 x 50 discs
Cefaclor 30	E113049	1 x 50 discs
Cefadroxil 30	E113057	1 x 50 discs
Cefamandole 30	E113030	1 x 50 discs
Cefepim 30	E113002	1 x 50 discs
Cefepim/Clavulan Acid 30/10	E113125	1 x 50 discs
Cefixime 5	E113058	1 x 50 discs
Cefoperazone 30	E113055	1 x 50 discs
Cefoperazone 75	E113074	1 x 50 discs
Cefoperazone/Sulbactam 3:1	E115727	1 x 50 discs
Cefotaxime 30	E113018	1 x 50 discs
Cefotaxime 5	E113064	1 x 50 discs
Cefotaxime/Clavulan Acid 30/10	E113092	1 x 50 discs
Cefoxitin 30	E113133	1 x 50 discs
Ceftazidime 10	E115629	1 x 50 discs
Ceftazidime 30	E113013	1 x 50 discs
Ceftazidime/Clavulan Acid 30/10	E113098	1 x 50 discs
Ceftibuten 30	E113072	1 x 50 discs
Ceftriaxone 30	E113056	1 x 50 discs
Cefuroxime 30	E113050	1 x 50 discs
Cephalexin 30	E113132	1 x 50 discs
Cephalothin 30	E113012	1 x 50 discs
Cephazolin 30	E113054	1 x 50 discs

Description	Cat.No	Package size
Cephadrine 30	E113066	1 x 50 discs
Chinopristina/Dalfopristina	E110104	1 x 50 discs
Chloramphenicol 30	E113016	1 x 50 discs
Ciprofloxacin 5	E113142	1 x 50 discs
Clarithromycin 15	E113060	1 x 50 discs
Clindamycin 2	E113022	1 x 50 discs
Cloxacilin 5	E113032	1 x 50 discs
Colistin sulphate 10	E113114	1 x 50 discs
Colistin sulphate 25	E113113	1 x 50 discs
Colistin sulphate 50	E115667	1 x 50 discs
Cotrimoxazole 23,75/1,25	E113122	1 x 50 discs
Doxycycline 30	E111412	1 x 50 discs
Ertapenem 10	E110502	1 x 50 discs
Erythromycin 15	E110515	1 x 50 discs
Fosfomycin 50	E115183	1 x 50 discs
Furazolidon 100	E115122	1 x 50 discs
Furazolidon 50	E115023	1 x 50 discs
Fusidic Acid 10	E110617	1 x 50 discs
Gentamicin 10	E110712	1 x 50 discs
Gentamicin 120	E110718	1 x 50 discs
Gentamicin 30	E115072	1 x 50 discs
Imipenem 10	E110913	1 x 50 discs
Kanamycin 30	E111112	1 x 50 discs
Levofloxacin 5	E110106	1 x 50 discs
Lincomycin 15	E115028	1 x 50 discs
Lincomycin 2	E111212	1 x 50 discs
Linezolid 10	E111246	1 x 50 discs
Linezolid 30	E110107	1 x 50 discs
Meropenem 10	E111352	1 x 50 discs
Metronidazole 5	E111328	1 x 50 discs
Moxifloxacin 5	E111364	1 x 50 discs
Mupirocin 5	E111330	1 x 50 discs
Nalidixic Acid 30	E111426	1 x 50 discs
Neomycin 10	E111436	1 x 50 discs
Neomycin 30	E111429	1 x 50 discs
Netilmicin 10	E111438	1 x 50 discs
Netilmicin 30	E111434	1 x 50 discs
Nitrofurantoin 100	E115034	1 x 50 discs

Antibiotic paper discs

Disques antibiogramme

Description	Cat.No	Package size
Nitrofurantoin 300	E111432	1 x 50 discs
Norfloxacin 10	E111433	1 x 50 discs
Novobiocin 30	E111465	1 x 50 discs
Novobiocin 5	E111463	1 x 50 discs
Ofloxacin 5	E111516	1 x 50 discs
Oxacillin 1	E111515	1 x 50 discs
Oxacillin 5	E111532	1 x 50 discs
Oxytetracyclin 30	E111522	1 x 50 discs
Pefloxacin 5	E111634	1 x 50 discs
Penicillin G 1 units	E111606	1 x 50 discs
Penicillin G 10 units	E111614	1 x 50 discs
Phenoxymethylpenicillin 10	E191712	1 x 50 discs
Piperacillin 100	E111616	1 x 50 discs
Piperacillin 30	E111618	1 x 50 discs
Piperacillin/Tazobactam 100/10	E110103	1 x 50 discs
Piperacillin/Tazobactam 75/10	E112036	1 x 50 discs
Piperacillin/Tazobactam 30/6	E191612	1 x 50 discs
Polymyxin B 300 units	E111632	1 x 50 discs
Rifampicin 5	E111816	1 x 50 discs

Description	Cat.No	Package size
Spectinomycin 100	E111922	1 x 50 discs
Spiramycin 100	E110113	1 x 50 discs
Streptomycin 10	E111925	1 x 50 discs
Streptomycin 25	E111926	1 x 50 discs
Streptomycin 300	E111946	1 x 50 discs
Sulphamethoxazole/Trimeth. 23,75/1,25 s.Cotrimoxazol	E113122	1 x 50 discs
Teicoplanin 30	E112022	1 x 50 discs
Telithromycin 15	E112054	1 x 50 discs
Tetracycline 30	E112010	1 x 50 discs
Ticarcillin 75	E112031	1 x 50 discs
Ticarcillin/Cavulan Acid 75/10	E112033	1 x 50 discs
Tigecycline 15	E112056	1 x 50 discs
Tobramycin 10	E112020	1 x 50 discs
Trimethoprim 1,25	E112042	1 x 50 discs
Trimethoprim 5	E112026	1 x 50 discs
Vancomycin 30	E112023	1 x 50 discs
Vancomycin 5	E112028	1 x 50 discs
Voriconazole 1	E100634	1 x 50 discs

Antifungal agents

Description	Cat.No	Package size
Amphotericin B 20	E110191	1 x 50 discs
Clotrimazol 10	E110391	1 x 50 discs
Econazol 10	E190512	1 x 50 discs
Fluconazol 10	E190622	1 x 50 discs
Fluconazol 25	E190632	1 x 50 discs
Ketoconazol 10	E191112	1 x 50 discs
Miconazol 10	E191312	1 x 50 discs
Nystatin 100 units	E191422	1 x 50 discs

Veterinary agents (not IVD)

Description	Cat.No	Package size
Cefacetril 30	E191816	1 x 50 discs
Cefquinome 30	E115729	1 x 50 discs
Enrofloxacin 5	E110510	1 x 50 discs
Florfenicol 30	E110619	1 x 50 discs
Flumequina 30	E115666	1 x 50 discs
Marbofloxacin 5	E100635	1 x 50 discs
Tylosin 30	E112015	1 x 50 discs

Equipment

Description	Cat.No	Package size
Discs Dispenser, 6 cartridges	EM0006	1 pcs

Identification

D'identification

Diagnostic discs

Disques d'identification

Description	Cat.No	Package size
Bacitracin 0,1	E171011	1 x 50 discs
Bacitracin 0,04 for Streptococci Group A identification	E171010	1 x 50 discs
BV Factor for Haemophilus spp. identification	E171042	1 x 50 discs
BVX Factor for Haemophilus spp. identification	E171052	1 x 50 discs
BX Factor for Haemophilus spp. identification	E171032	1 x 50 discs
Furazolidon 50 for differentiation Staphylococcus and Micrococcus	E115023	1 x 50 discs
Blanc discs	E110100	1 x 50 discs
Metronidazole 50 µg for G. vaginalis identification	E171061	1 x 50 discs
Novobiocin 5 for Staphylococcus saprophyticus identification	E111463	1 x 50 discs
Oxidase strips	10003324	1 x 50 strips
Optochin 5 µg for S. pneumoniae identification	E171021	1 x 50 discs
Trimethoprim 5 for G. vaginalis identification	E112026	1 x 50 discs
V + X Factor for Haemophilus spp. identification	E171051	1 x 50 discs
V Factor for Haemophilus spp. identification	E171041	1 x 50 discs
X Factor for Haemophilus spp. identification	E171031	1 x 50 discs

Latex tests

Tests d'agglutination latex *

Description	Cat.No	Package size
Prolex™ Staph Kit	PL.080B	100 tests
Prolex™ Staph Kit	PL.081B	300 tests
Prolex™ Staph Xtra Kit	PL.1080	100 tests
Prolex™ Staph Xtra Kit	PL.1081	300 tests
Prolex™ Strep Grouping Kit	PL.030	60 tests
Prolex™ E.coli O157 Kit	PL.070B	50 tests
Prolex™ E.coli O157 Kit	PL.071B	100 tests
Prolex™ E. coli non O157 Kit	PL.1070	50 tests
Prolex™ L.pneumo. SG1 DFA Kit	PL.241	50 tests
Prolex™ L.pneumo. SG1-14 DFA Kit	PL.242	50 tests

*don't shipp to all countries

Reagents and tests

Reactifs et tests*

Description	Cat.No	Package size
CANDIFAST® - Test combining identification and antifungal resistance testing of the main pathogenic yeasts	44008	8 tests
CANDIFAST® - Test combining identification and antifungal resistance testing of the main pathogenic yeasts	44030	30 tests
CANDIFAST® ES TWIN - Antifungal resistance testing of yeasts responsible for superficial mycosis	44130	30 tests
FUNGIFAST® AFG - Antifungal susceptibility testing of the main pathogenic yeasts	44412	12 tests
FUNGIFAST® - Test combining identification and antifungal resistance testing of the main pathogenic yeasts	44408	8 tests
FUNGIFAST® - Test combining identification and antifungal resistance testing of the main pathogenic yeasts	44430	30 tests
MYCOFAST® RevolutionN - Detection, enumeration, identification and complete antimicrobial susceptibility testing of urogenital mycoplasma: DIRECT METHOD	00060	25 tests
Mycoplasma Stabilizer	00064	2 x 5 ml
Bile salts	464/T	1 x 5 ml
Coagulase Plasma	PL.850-3	10 x 3 ml
EDTA 0,5 M	462/T	1 x 5 ml
Indole Reagent	PL.391-10	10 ml
Kovacs Reagent	PL.375	100 ml
McFarland Standard - Set of 0.5-4.0	SD 2350	5 x 10 ml
McFarland Standard 0.5	SD 2300	10 ml
McFarland Standard 1.0	SD 2301	10 ml
McFarland Standard 2.0	SD 2302	10 ml
McFarland Standard 3.0	SD 2303	10 ml
McFarland Standard 4.0	SD 2304	10 ml
Oxidase Reagent	PL.390	15 ml
Phenylboronic acid	461/T	1 x 5 ml

*don't shipp to all countries

Stains

Colorains

Description	Cat.No	Package size
Gram Staining Kit	PL.8055/25	4 x 250 ml
Gram Staining Kit	PL.8056/25	4 x 250 ml
Gram Staining Kit	PL.8057/25	4 x 250 ml
TB Staining Kit	PL.8060/25	4 x 250 ml
TB Staining Kit	PL.8061/25	4 x 250 ml
Crypto. Staining Kit	PL.8060	7 x 250 ml
Immersion Oil	PL.396	50 ml

Cryogenic storage of microorganisms – VIABANK

stockage cryogénique des microorganismes - VIABANK

VIABANK™ vials containing a minimum of 25 glass beads covered in a cryopreservative solution for long term storage of microorganisms at low temperatures. Individual beads can be removed while frozen at any time for establishing fresh cultures. Each box contains 80 VIABANK vials with either 4 different coloured caps in each box or a choice of one colour from the five options.

Description	Cat.No	Package size
VIABANK - Vials with red beads	VIRA	80 vials
VIABANK - Vials with green beads	VIGA	80 vials
VIABANK - Vials with yellow beads	VIYA	80 vials
VIABANK - Vials with blue beads	VIBA	80 vials
VIABANK - Vials with multi coloured beads	VIMA	80 vials

Cell culture media and buffers

Millieux de culture de cellules et

Description	Cat.No	Package size
Buffered Salt Dilution (PBS)	BT 5105.01	1 x 100 ml
Buffered Salt Dilution (PBS)	BT 5105.02	1 x 200 ml
Buffered Salt Dilution (PBS)	BT 5105.05	1 x 500 ml
PBS with calcium and magnesium ions	BT 5290.01	1 x 100 ml
PBS with calcium and magnesium ions	BT 5290.02	1 x 200 ml
PBS with calcium and magnesium ions	BT 5290.05	1 x 500 ml
RPMI 1640 + NaHCO ₃ + L-Glutamina + phenol red	BT 5295.01	1 x 100 ml
RPMI 1640 + NaHCO ₃ + L-Glutamina + phenol red	BT 5295.02	1 x 200 ml
RPMI 1640 + NaHCO ₃ + L-Glutamina + phenol red	BT 5295.05	1 x 500 ml



BioMaxima S.A.
ul. Vetterów 5
20-277 Lublin, Poland
www.biomaxima.com

phone +48 81 467 90 61
fax +48 81 744 29 15
export@biomaxima.com
support@biomaxima.com
microbiology@biomaxima.com